




generalmedical



2024

**Advanced Dental, Surgical and
Implant Product Specialists**



PIEZOSURGERY®

PIEZOSURGERY® TOUCH & WHITE EXPERIENCE PERFORMANCE



FINANCE


PRICE
DEALS



Extraction Kit



Molar Extraction Kit




Explantation Kit



Osteotomy Kit



Resective Perio Kit



Retro Surgical Kit



Piezo Crestal Lift Kit



Sinus Lift Kit



Sinus Lift Lateral Kit

Order Online at www.generalmedical.co.uk

Or call 01380 734990 or email info@generalmedical.co.uk




PIEZOSURGERY & MOTORS

6



BONE MANAGEMENT

10




INSTRUMENTS

20



BIOMATERIALS

42



IMPLANTS

64

LOCATORS AND
RESTORATIVE OPTIONS

68



STERILE SUPPLIES

76

Choukroun Quattro PRF System

A new concept in Smart Blood Generation, reduces the need for expensive xenografts and membranes.

SEE PAGE 52



BiLumix
SHADOWLESS HEADLAMP

Illuminate Your Vision

For any procedure to illuminate an area via a unique, dual light source that eliminates any shadow.



BiLumix

SHADOWLESS HEADLAMP

Headlamps Features


One-touch On-Off

Additional on and off switch on back to help you control your light more efficiently during your procedures.



Power Indicator

Easy way to check battery status.



Easy Control

Quickly change brightness of LED with dial control.



Loupes Specifications



1.5x:

Weight: 18.8g

Various Working Distance



2.5x:

Weight: 48.7g

4.0x:

Weight: 67.9g

Adjustable Working Distance
(25 – 60cm)

Piezosurgery Touch Cutting edge technology



**From the world leaders
in Piezosurgery®**

**All insert tips purchased
with new units will be
discounted by 15%**

Piezosurgery White

All of the benefits - at a lower price!



**Call now for
latest offers**



MEC-03540011
Trolley with three
toughened glass
shelves and power
outlet.

Piezosurgery Insert Tip Kits


All Mectron Piezosurgery insert tips are available as separate items.
Please contact us for more details.



MEC-01520003
Extraction Kit



MEC-01520025
Molar Extraction Kit




MEC-01520021
Explantation Kit



MEC-01520010
Osteotomy Kit



MEC-01520022
Resective Perio Kit



MEC-01520005
Retro Surgical Kit



MEC-01520024
Piezo Crestal Lift Kit



MEC-01520002
Sinus Lift Kit




MEC-01520023
Sinus Lift Lateral Kit


Piezosurgery Insert Tips

Our range at a glance


Osteotomy




Osteotomy








Explantation



Osteoplasty



*only available for PIEZOSURGERY®3, PIEZOSURGERY® touch and PIEZOSURGERY® white

Implant Site Preparation**Extraction****Third Molar Extraction****Membrane Detachment****Periosteal Elevation****Endodontics****Periodontology****Implant Cleaning**

Membrane Fixation Set

Tack system starter set

Tack holder, insertion instrument, removal instrument -
10 x 3mm tacks
ZEP-TACKSET

For secure positioning of all types of membranes including Titanium mesh

ZEP-47.520.01
Pin Applicator

TN-002 Pack of 10
3mm Titanium Tacks

ZEP-47.847.12
Pin Remover

TN-003 Pack of 10
5mm Titanium Tacks

ZEP-85.26.00
Tack Holder



Krekeler Sliding Caliper

These sliding calipers combine a variety of functions in one instrument, facilitating precise positioning of implants. By popular request Zepf has made a modification to the existing Krekeler sliding caliper - a locking screw to fix the space measured.

ZEP-31.693.10



ZEP-47.958.00
Bone Mill Forceps



Witzel All Purpose Forceps

Diamond tipped root splinter forceps for upper
or lower jaw.

ZEP-14.700.01



HELMUT ZEPP
MEDIZINTECHNIK GMBH

Master-Mill

Surgical Bone Mill

The Master-Mill is a surgical bone mill used for crushing autologous bone blocks and thus preserving small autologous bone chips. With the help of these bone chips, bone defects of all shapes can be optimally filled. The individual components of the Master-Mill are made of high quality, surgical stainless steel. All parts feature an optimal fit and a stable geometry, allowing the Master-Mill to be easy to handle with ease of assembly and intuitive use.

MEI-BKM00



NEW

Mills, Clamps and Files

Bone Block Holder

ZEP-23.105.50



Mini Bone Mill

ZEP-47.955.00



Beyer Rongeur

Double action bone rongeur

ZEP-42.650.18



Master-Pin-Control

Pin System for Membrane Fixation

Developed with Dr. Istvan Urban

Master-Pin-Control is designed especially for the fixation of resorbable, non-resorbable and titanium reinforced membranes. The included pins consist of titanium alloy Grade 5 and have a groove which increases the surface of the pins. Due to the groove and a unique head design, the pins can be loosened easily and safely which is especially important after a successful healing period. With their sharp tip and very stable shank, they can easily and precisely be inserted into dense cortical bone.

35x Pins

MEI-BMP00



	Initial burs 2 pcs.	Twist drills each 2 pcs.	Pins 35 pcs.	
Fig.	186RF	203RF	203RF	MP10*
Shank ¹	204	205	205	-
Size ²	018	006	008	-
Length mm	12.0	7.0	7.0	3.65
	-	-	-	0.87
	1.8	0.6	0.8	0.95
Head-Ø mm	-	-	-	2.51

¹ 204=RA, 205=RA L

² Largest working part diameter in 1/10 mm

Minimal diameter External diameter



Optionally available:

Twist Drill with Stop
for decorticating

203S



Refill Pins:

Pins 10 (MP10*)

Pins 50 (MP50*)

Pins 100 (MP100*)



Screwdriver for Pin




Fig.

MP11

illustrated 1:1 (Length mm: 35.0)

Pin holder



Fig.

MP12

illustrated 1:2 (Length mm: 170.0)

Fixing holder



Fig.

MP14

illustrated 1:2 (Length mm: 165.0)

Micro Screw System

Olsberger Konzept by Prof. Dr. Fouad Khoury



The Bone Management® Micro Screw Systems are used for the safe fixation and stabilisation of cortical bone grafts. The special feature of these systems are the osteosynthesis screws, which have a diameter of 1.0 mm and 1.2 mm. These Micro Screws are made of surgical stainless steel, which gives them high stability despite their reduced diameter.

MEI-BMS00

Optionally available:

915DM medium NEW

Fig.	Shank	Ref.-No.	L mm
915DM	HP	806 104 35M 524	220
	unmounted	806 900 35M 524	220

Super flexible

















	Pilot burs		Initial bur		Screw driver	Claw for manual screwdriver	Screw driver RA
Fig.	MSPB1	MSPB2	202RF	202RF	MSSDM	MSGSD	MSSDW
Shank ¹	204	204	206	206	-	-	204
Size ²	008	010	008	010	-	-	-
Length mm	14.0	14.0	14.0	14.0	142.2	33.0	24.0
Minimal diameter	-	-	-	-	-	-	-
External diameter	0.8	1.0	0.8	1.0	-	-	-

¹ 204 = RA, 206 = RA XL

² Largest working part diameter in 1/10 mm

Minimal diameter External diameter

	Screw (2 pcs.)	Screws (each 3 pcs.)			Screws (each 2 pcs.)		Screws (1 pcs.)	Screws (each 2 pcs.)			Screws (each 1 pcs.)	
												
Fig.	MSS10 040*	MSS10 060*	MSS10 080*	MSS10 100*	MSS10 120*	MSS10 140*	MSS12 040*	MSS12 060*	MSS12 080*	MSS12 100*	MSS12 120*	MSS12 140*
Shank ¹	-	-	-	-	-	-	-	-	-	-	-	-
Size ²	-	-	-	-	-	-	-	-	-	-	-	-
Length mm	4.0	6.0	8.0	10.0	12.0	14.0	4.0	6.0	8.0	10.0	12.0	14.0
	0.6	0.6	0.6	0.6	0.6	0.6	0.7	0.7	0.7	0.7	0.7	0.7
	1.0	1.0	1.0	1.0	1.0	1.0	1.2	1.2	1.2	1.2	1.2	1.2
Head-Ø mm	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7	1.7

² Largest working part diameter in 1/10 mm

Minimal diameter

External diameter

* (Symbol)

Alveoplasty Surgical Kit










Surgical Contouring and Modification Kit





The Alveoplasty Surgical Kit contains the most important instruments for surgical contouring and modification of the alveolar ridge when unwanted unevenness of the jawbone occurs after tooth extraction or bone degeneration. In this case, the kit offers the best way to smooth the jawbone. The Alveoplasty Kit contains two initial burs and two Lindemann burs, which allow for the performance of precise, initial osteotomies, even in hard cortical bone. To perform the alveoplasty itself, five different carbide cutters and a round diamond are included. The round working parts optimally protect the soft tissue against injuries. Thus, the alveolar ridge is prepared carefully and precisely for implant insertion or denture estoration.



MEI-ALV18

	Initial burs (2 pcs.)			Surgical Cutters		Tungsten Carbide Instruments				Diamond Instrument
										
Fig.	I86RF	HM166**	I65RF	HM71**	HM71**	HM72GX**	HM77GX**	HM251GX**	801*	
Shank ¹	204	205	205	104	104	104	104	104	204	
Size ²	018	021	023	040	050	060	060	040	033	
Length mm	12.0	11.0	7.0	-	-	12.0	11.5	8.3	-	
* -	1.8	2.1	2.3	4.0	5.0	6.0	6.0	4.0	3.3	

¹ 104=HP, 204=RA, 205=RA L ² Largest working part diameter in 1/10 mm * External diameter *  ** 

Diamond Burs



For surgical (straight) handpiece:

MEI-801-104-021	Medium grit (2.1mm)
MEI-801-104-023	Medium grit (2.3mm)
MEI-801-104-027	Medium grit (2.7mm)
MEI-801-104-033	Medium grit (3.3mm)
MEI-801-104-050	Medium grit (5mm)
MEI-801G-104-021	Coarse grit (2.1mm)
MEI-801G-104-023	Coarse grit (2.3mm)
MEI-801G-104-033	Coarse grit (3.3mm)
MEI-801G-104-050	Coarse grit (5mm)

For contra-angle (latch grip) handpieces:

MEI-801-204-021	Medium grit (2.1mm)
MEI-801-204-023	Medium grit (2.3mm)
MEI-801-204-027	Medium grit (2.7mm)
MEI-801-204-033	Medium grit (3.3mm)

Master-Core Professional

Developed by Prof. Istvan Urban

The Master-Core Professional is used for the simple and safe extraction of bone cylinders in the context of bone augmentation. For this purpose, the system contains trephines with different lengths and diameters. Due to the different lengths, the trephines have an automatic depth stop, which ensures the protection of anatomical structures and thus allows for a particularly safe and controlled bone extraction. In addition, depth markings on the trephines' working parts enable additional depth control during the extraction. The black coating of the trephines allows for glare-free work and the optimal visibility of the depth markings. To be able to respond more precisely in specific indications, the Master-Core Professional contains a broad range of instruments.

MEI-BMCPR



	Initial bur Trephines				Drill Extension	
Fig.	DV229	DV229	DV229	DV229	BV024	BV025
Shank ¹	204	204	204	204	205	205
Size ²	040	050	060	070	-	-
Length mm	2.0	2.0	2.0	2.0	23.5	25.5
Internal diameter	4.0	5.0	6.0	7.0	-	-
External diameter	4.7	5.7	6.7	7.7	-	-

¹ 204=RA, 205=RA L

² Largest working part diameter in 1/10 mm

Internal diameter

External diameter

	Trephines									
Fig.	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU
Shank ¹	206	206	206	206	206	206	206	206	206	206
Size ²	503	504	505	506	507	603	604	605	606	607
Length mm	3.0	4.0	5.0	6.0	7.0	3.0	4.0	5.0	6.0	7.0
Internal diameter	4.0	4.0	4.0	4.0	4.0	5.0	5.0	5.0	5.0	5.0
External diameter	5.0	5.0	5.0	5.0	5.0	6.0	6.0	6.0	6.0	6.0

¹ 206=RA XL

² Largest working part diameter in 1/10 mm

Internal diameter

External diameter

	Trephines									
Fig.	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU	229IU
Shank ¹	206	206	206	206	206	206	206	206	206	206
Size ²	703	704	705	706	707	803	804	805	806	807
Length mm	3.0	4.0	5.0	6.0	7.0	3.0	4.0	5.0	6.0	7.0
Internal diameter	6.0	6.0	6.0	6.0	6.0	7.0	7.0	7.0	7.0	7.0
External diameter	7.0	7.0	7.0	7.0	7.0	8.0	8.0	8.0	8.0	8.0

¹ 206=RA XL

² Largest working part diameter in 1/10 mm

Internal diameter

External diameter


Disposable tissue punches

SAL-015 3.5mm box of 20
SAL-015-1 3.5mm box of 10
SAL-016 4.0mm box of 20
SAL-016-1 4.0mm box of 10
SAL-017 5.0mm box of 20
SAL-017-1 5.0mm box of 10
SAL-039 6.0mm box of 20
SAL-039-1 6.0mm box of 10



Easy Removal of Soft Tissue

ZEP-08.921.00



231 DC extra fine

Fig.	Shank	Ref.-No.		1	1	1
			L mm	0,3	0,3	0,3
231DC	RA	806 204 064 504		070*	100*	130*



* im Saw Basic Kit 7150 enthalten
* contained in the Saw Basic Kit 7150
* figuran en el Saw Basic Kit 7150

Bone Chisel Set

ZEP-47.949.95 complete set
Includes all 6 chisels listed plus storage cassette.



Pointed



Flat



ZEP-47.949.11 (4mm)



ZEP-47.949.12 (6mm)



ZEP-47.949.13 (8mm)



ZEP-47.949.21 (4mm)



ZEP-47.949.22 (6mm)



ZEP-47.949.23 (8mm)



Mallets

ZEP-41.500.17
Lead filled mallet



ZEP-41.501.01
Lightweight handle, exchangeable
25mm teflon inserts



ZEP-41.503.01
Ferrozell mallet

Mallet made of a revolutionary wood-like plastic material giving softer more tactile handling characteristics. New 250 mm length fits all standard cassettes.



Membrane Block Fixation and Tenting Screws

All in one box!

The **Pro-fix™ Precision Fixation System** is manufactured to precise tolerances to ensure easy pick-up of screws, stable transfer to the surgical site, and quick engagement in cortical bone. All kit components are labelled, organised, stored and sterilised together to simplify things for the surgeon and surgical team.



Driver Handle and Blade

The stainless steel driver handle and cruciform driver blade fit together to make up a drive system that ensures easy pick-up of screws and stable transfer to the surgical site.

Tecapro Storage Tray

The autoclavable Tecapro storage tray is cleverly designed to store not only fixation instruments, but also a variety of membrane fixation, bone fixation and tenting screws needed in dental bone grafting.

Screw Organiser Dial

Up to 100 fixation screws - including bone fixation, membrane fixation, and tenting - may be stored in a built-in screw organiser dial designed and labelled for easy identification, easy storage, and simple re-ordering.

Membrane screw



Tenting screw



Bone screw



Pro-fix™ Membrane Fixation Screws are designed as an attractive alternative to using tacks for membrane stabilisation. Easy pick-up, solid stability of the screw during transfer to the surgical site, and easy placement make membrane fixation fast and easy.

Membrane Fixation Kit

CYT-PFMK20

- 1 Storage tray with screw organiser dial
- 1 Stainless steel driver handle
- 1 56mm cruciform driver blade
- 1 76mm cruciform driver blade
- 20 1.5mm x 3.0mm self-drilling titanium alloy membrane fixation screws.

**Ask about our Pro-fix Combi Kit which includes Membrane, Block and Tenting screws*

Self-Drilling Membrane Fixation Screws

CYT-PFMF-5 5 screws

Individual Components

Blades work universally with all Pro-Fix screws

CYT-PFDH

Stainless steel driver handle

CYT-PFDB

76mm cruciform driver blade

CYT-PFDB56

56mm cruciform driver blade

CYT-PFDBCA

Contra angle blade - (24mm long; 10mm exposed distal length)

CYT-PFPD

1.2mm diameter latch type pilot drill

Pro-fix™ Tenting Screws are designed with a self-drilling tip, polished neck, and broader head to maintain space under resorbable and non-resorbable membranes in horizontal and vertical bone regeneration procedures.

Tenting Kit

CYT-PFTK12

- 1 Storage tray with screw organiser dial
- 1 Stainless steel driver handle
- 1 56mm cruciform driver blade
- 1 76mm cruciform driver blade
- 4 1.5mm x 3.0mm self-drilling titanium alloy tenting screws. *(3.0mm polished neck + 4.0mm threaded portion = 7mm total length)*
- 4 1.5mm x 4.0mm self-drilling titanium alloy tenting screws. *(4.0mm polished neck + 4.0mm threaded portion = 8mm total length)*
- 4 1.5mm x 5.0mm self-drilling titanium alloy tenting screws. *(5.0mm polished neck + 4mm threaded portion = 9mm total length)*

Self-Drilling Tenting Screws

CYT-PFT3-5 5 screws

3.0mm polished neck + 4.0mm threaded portion = 7mm total length

CYT-PFT4-5 5 screws

4.0mm polished neck + 4.0mm threaded portion = 8mm total length

CYT-PFT5-5 5 screws

5.0mm polished neck + 4.0mm threaded portion = 9mm total length

CYT-PFT8 1 screw

8.0mm fully threaded tenting screw

CYT-PFT10 1 screw

10.0mm fully threaded tenting screw

Pro-fix™ Bone Fixation Screws

are designed with finer pitched, self-tapping threads that give the screws greater clamping force while using less driver torque. The screws' threads are equipped with a cutting flute that allows for easier insertion into harder bone. The screws are placed into a 1.2mm pre-drilled pilot hole.

Bone Fixation Kit

CYT-PFBK12

- 1 Storage tray with screw organiser dial
- 1 Stainless steel driver handle
- 1 56mm cruciform driver blade
- 1 76mm cruciform driver blade
- 1 1.2mm diameter latch type pilot drill
- 2 1.5mm x 8mm bone fixation screws
- 4 1.5mm x 10mm bone fixation screws
- 2 1.5mm x 14mm bone fixation screws

MEI-B2001

Lagdrill 1.6mm diameter ensures screw passes easily through the block.

Self-Tapping Bone Fixation Screws

CYT-PFB8-5 5 screws

1.5mm x 8mm

CYT-PFB10-5 5 screws

1.5mm x 10mm

CYT-PFB12-5 5 screws

1.5mm x 12mm

CYT-PFB14-5 5 screws

1.5mm x 14mm

MEI-B2001

Lagdrill 1.6mm diameter ensures screw passes easily through the block.



ZEP-12.300.08

Benex[®] II Extractor
incl. Support Disc



ZEP-12.300.80

Quadrant Support
to bridgeover bigger gaps
and for universal molding

Benex[®] II The individual components

acc. to Dr. med., med. dent. Benno Syfrig

ZEP-12.302.00 Benex[®] II Extraction System - with full size cassette to add other instrument to

ZEP-12.303.00 Benex[®] II Basic Kit - with ½ size cassette

Apart from the Washbasket 1/1 with Lid and Press Button Lock (**85.195.10**) incl. the Benex II Tray / Rack (**12.302.01**) for the Extraction System and the Washbasket 1/2 with Lid and Press Button Lock (**85.194.10**) for the Basic Kit, the following components are included in both Sets (pict. scale 1:1):



ZEP-12.300.15

Replacement Support
Disc, 8 mm (PTFE)



ZEP-12.300.20

Pullrope, 48 mm,
PU: 2 pieces



ZEP-12.300.47

Driver Guide, short



ZEP-12.300.60

1 Screw 10 mm
Ø 1.6 mm, S = Short



ZEP-12.300.65

1 Screw 10 mm
Ø 2.1 mm, SF = Short & Fat



ZEP-12.300.30

1 Diamond coated Drill
for Screws Ø 1.6 mm
12.300.60 and 12.300.70



ZEP-12.300.70

1 Screw 16 mm
Ø 1.6 mm, L = Long



ZEP-12.300.75

1 Screw 16 mm
Ø 2.1 mm, LF = Long & Fat



ZEP-12.300.35

1 Diamond coated Drill
for Screws Ø 2.1 mm
12.300.65 and 12.300.75

OPTIONALLY AVAILABLE




ZEP-12.300.17

Support Disc,
diagonally left



ZEP-12.300.16

Support Disc,
diagonally right



ZEP-12.300.48

Driver Guide, long



ZEP-12.300.45

Blade for

ZEP-47.525.51

Driver Guide FD

Periotomes

For separation of Sharpey's fibres in the sulcus.

ZEP-003 Periotome set (includes handle and blades P1, P2 & P3)



ZEP-004 Blade P1



ZEP-005 Blade P2



ZEP-007 Periotome Handle



ZEP-006 Blade P3



ZEP-002 Blade P4

Power Periotomes

The Power Periotomes allow gentle loosening of the fibres in the sulcus.
The handle provides optimal power transmission and controlled luxation.



ZEP-26.690.01 Power Periotome P1 XL, purple



ZEP-26.690.02 Power Periotome P2 XL, mesial - blue



ZEP-26.690.03 Power Periotome P3 XL, distal - green

FLEX-EX Power Periotome

A symbiosis of Power Periotomes, Elevators and Xtool.


The flexible working tip will enable you to build up a phenomenal pressure for luxation in a radial direction and a perfect match to the contour of the tooth.



ZEP-26.690.10 FLEX-EX Power Periotome

Xtool Luxators

Top quality 'Hildebrand' Luxators.
All items available separately.



ZEP-17.007.00

Complete Set contains
5 luxators:

ZEP-17.007.01*

ZEP-17.007.02*

ZEP-17.007.03*

ZEP-17.007.04*

ZEP-17.007.05*

plus ZEP-41.855.01Z curette plus
free round Arkansas Stone.



ZEP-17.007.01*
ZEP-17.007.01TI



ZEP-17.007.02*
ZEP-17.007.02TI



ZEP-17.007.03*
ZEP-17.007.03TI



ZEP-17.007.04*
ZEP-17.007.04TI



ZEP-17.007.05*
ZEP-17.007.05TI



ZEP-17.007.06
ZEP-17.007.06TI



ZEP-17.007.07
ZEP-17.007.07TI



ZEP-17.008.01
ZEP-17.008.01TI



ZEP-17.008.02
ZEP-17.008.02TI



ZEP-17.008.03
ZEP-17.008.03TI



ZEP-24.950.00 Stone Oil

ZEP-24.923.10 Arkansas Stone (round) for sharpening



ZEP-17.008.04
ZEP-17.008.04TI



ZEP-17.008.05
ZEP-17.008.05TI



ZEP-17.008.06
ZEP-17.008.06TI



ZEP-17.008.08
ZEP-17.008.08TI



ZEP-17.008.09
ZEP-17.008.09TI



ZEP-17.008.10
ZEP-17.008.10TI



ZEP-17.008.11
ZEP-17.008.11TI



ZEP-17.008.90
ZEP-17.008.90TI



All available in
Onyx (black),
TI codes

RoBa Extraction Forceps

RoBa Range

RoBa deep grip forceps have concave beaks developed to ensure maximal adaptation on convex shaped teeth. Due to their special shape, these forceps assure parallel and maximum grip of the crown in any situation. Root fracture is avoided. Also, no tilting of the tooth during rotational or lifting movements prevents damage to the buccal plate.

Also available in diamond tipped for better grip and Onyx (black) non-reflective high-grip coating.



Diamond Version

All extraction forceps are available with diamond coating for a better grip. Indicated by gold plating on the handle.



ZEP-12.234.07Z
ZEP-12.234.07ZD
ZEP-12.234.07ZTI

Stainless Steel
Diamond coated
Onyx coated

Upper 3,2,1
Narrow



ZEP-12.234.08Z
ZEP-12.234.08ZD
ZEP-12.234.08ZTI

Upper 3,2,1
Wide



ZEP-12.235.07Z
ZEP-12.235.07ZD
ZEP-12.235.07ZTI

Upper 5,4
Narrow



ZEP-12.235.08Z
ZEP-12.235.08ZD
ZEP-12.235.08ZTI

Upper 5,4
Wide



ZEP-12.217.00Z
ZEP-12.217.00ZD
ZEP-12.217.00ZTI

Upper Right
6,7



ZEP-12.218.00Z
ZEP-12.218.00ZD
ZEP-12.218.00ZTI

Upper Left
6,7



ZEP-12.267.01Z
ZEP-12.267.01ZD
ZEP-12.267.01ZTI

Upper 8

ZEPF
roba
 EDITION
 PATENTED

HELMUT ZEPF
 MEDIZINTECHNIK GMBH

Titanium Aluminium Nitride Coating

Users with the highest demands will appreciate the elegant, black finish providing the instrument with a non-reflecting, extremely smooth and scratch-resistant surface.

onyx
 DESIGN



ZEP-12.279.90Z
ZEP-12.279.90ZD
ZEP-12.279.90ZTI

Lower 8



ZEP-12.236.07Z
ZEP-12.236.07ZD
ZEP-12.236.07ZTI

Lower 3,2,1
Narrow



ZEP-12.236.08Z
ZEP-12.236.08ZD
ZEP-12.236.08ZTI

Lower 3,2,1
Wide



ZEP-12.213.00Z
ZEP-12.213.00ZD
ZEP-12.213.00ZTI

Lower 5,4



ZEP-12.222.00Z
ZEP-12.222.00ZD
ZEP-12.222.00ZTI

Lower 7,6
90° Pointed



ZEP-12.221.90Z
ZEP-12.221.90ZD
ZEP-12.221.90ZTI

Lower 7,6
offset, Flat end
Pointed

Periosteal Elevators

Buser Type Elevators

Each instrument has a pointed end and a rounded working end. The pointed working end is designed to lift the papilla before the periosteum is elevated with the rounded working end.



ZEP-41.864.30 Rounded end 2.8mm lance 2.4mm



ZEP-41.864.40 Rounded end 4.2mm lance 3.8mm



ZEP-41.864.50 Rounded end 4.2mm lance 3.8mm (Slight bend on lance and curette ends)



ZEP-41.864.13 Zepf line micro elevator 2.6mm ends.

Molt Elevators



No. 1 ZEP-41.862.11



No 2 ZEP-41.862.12



No 3 ZEP-41.862.13



No 4 ZEP-41.862.20

Prichard



ZEP-41.878.11

Seldin



ZEP-41.877.23 Backbone of periosteal elevators, excellent for general flap retraction.

Papillex Elevator

For elevations around the papilla, like a Buser but finer, exchangeable tips.
ZEP-46.035.05

All Purpose Elevator

ZEP-24.747.03

Allen Elevator

ZEP-24.437.11

Papilla Elevator

ZEP-46.035.00 General purpose ultrafine elevator.



ZEP-24.747.05 Curette - Double ended, ultra-fine

Lucas Currettes

ZEP-41.855.00Z - 1.5mm double ended

ZEP-41.855.01Z - 2.5mm double ended

ZEP-41.855.02Z - 3.0mm double ended

ZEP-41.855.04Z - 5.0mm double ended

For serrated end version of these currettes add suffix 'S' to the code.

Freer

ZEP-41.864.02

Mitchells Trimmer

ZEP-24.353.02

Supercut Operating Scissors

Zepf manufacture a huge range of operating scissors including titanium nitride coated long life versions. Please call for the latest Zepf catalogue.

Iris Scissors

Extra large rings, 115mm length

ZEP-46.050.11SC Straight

ZEP-46.051.11SC Curved




Goldman-Fox

130mm, tapered point

ZEP-46.221.12

ZEP-46.200.13SC Straight

ZEP-46.201.13SC Curved



Joseph Scissors

160mm, curved

ZEP-46.081.16SC

ZEP-46.081.16TIC



Le Grange

115mm double curved for easier access to posterior areas

ZEP-46.057.11SC



Suture Removal Scissors

Beck curved micro scissors with blade notch to cut and pull sutures. Fantastic product for all fine control tissue procedures.

ZEP-46.319.17 Micro serrated suture removal scissors

ZEP-46.319.17N Micro serrated suture removal scissors with notch





Onyx Scissors

The Onyx coating offers 3-5 times higher surface hardness compared to traditional scissors. In combination with 'Supercut' grinding, this guarantees an extremely long product life as well as very high precision.

ZEP-46.051.11TISC
Iris scissors - curved

ZEP-46.050.11TISC
Iris scissors - straight



ZEP-46.057.11TISC
La Grange scissors - double curved



ZEP-46.081.16TISC
Joseph scissors - curved and serrated



ZEP-46.201.13TISC
Goldman Fox scissors - curved



ZEP-46.319.17
Microsurgery scissors - curved 17cm

Also available

ZEP-46.319.17TISC
Onyx Microsurgery scissors -
curved 17cm with extra sharp long life nitride coating



Needle Holders

Castroviejo

Micro needle holders (5/0 - 8/0) with Spin-Lock hinge technology to prevent suture snags.

ZEP-41.010.17TC Straight

ZEP-41.011.17TC Curved

ZEP-41.010.17D Diamond Straight

ZEP-41.011.17D Diamond Curved



Macro Needle Holder

4/0 to 6/0

ZEP-41.010.17TC-M



Crile Wood

150mm (3/0 - 4/0)

ZEP-41.252.15TC

also available in 180mm

ZEP-41.252.18TC




Swedish Pattern

150mm (4/0 - 5/0)


A fine needle holder with flexible shafts to increase tactile feel and dampen locking mechanism.

ZEP-41.246.15TC



Dr Weng Hold 'n Cut Micro Needle Holders

ZEP-41.015.17



Tweezers

De Bakey

160mm
ZEP-50.226.16

College Tweezers

150mm
ZEP-22.011.03

MicroAdson

Toothed 150mm
ZEP-22.488.15

MicroAdson

Serrated 150mm, 1mm tips
ZEP-22.272.15

MicroAdson

Modified 150mm, toothed,
0.6mm fine tips
ZEP-22.489.01

Beck Atraumatic Forceps

Special shaped working tips allow
for gentle grip and minimal trauma
to tissues - a fantastic all purpose
forcep
ZEP-22.824.17

Micro Suture Forceps

ZEP-22.814.17

Atraumatic Tissue Forceps

ZEP-22.812.17

Cooley micro forceps -
serrated

ZEP-22.815.17 (straight)

ZEP-22.816.17 (curved)

ZEP-22.823.17 Extended reach
micro forceps, serrated, for pos-
terior region

ZEP-22.810.17

Micro forceps 17.5cm straight
(with teeth)



Retractors



Cawood-Minnesota

This is the go-to retractor for all surgical and implant work, an essential for every surgical kit.

ZEP-37.437.15



Cawood-Minnesota Modified

For longer procedures this modified Cawood is ideal, smoother edges make it easier to grip and more comfortable to hold.

ZEP-37.437.17



Kim Retractor

Designed to retract soft tissue flaps in the posterior region particularly during wisdom tooth removal. The working end is angled at 45 degrees, rounded and serrated so doesn't slip.

ZEP-37.446.04



Weider Retractor

Cheek retractor and tongue depressor

ZEP-37.426.01 small

ZEP-37.426.02 large



Bowdler-Henry's

Rake retractor for impacted lower wisdom teeth

ZEP-38.452.00



Laster Retractor

For impacted upper wisdom teeth, retract cheek, mucoperiosteum and hold an instrument all with the same hand.


ZEP-37.438.15



Vestibulum Retractor

Very useful in large augmentations upper jaw and lateral sinus procedures.

ZEP-38.448.50



Maty

Cheek and lip retractor

For upper and lower jaws, lightweight, no trauma - an excellent all round retractor

ZEP-37.454.00 adult

ZEP-37.454.02 paediatric

Scalpels and Blades



Miniature Blade Handle

Designed to make blade fitting and removal easier.
SM-SF13



Zepf Miniature Blade Handle

Ultra lightweight and counter balanced.
ZEP-46.013.00Z



Lightweight Round Handle

Adaptable handle for extra control, standard blades.
SM-B3



Standard Flat Handle

ZEP-46.010.03



Zepf Round Handle

Standard blades, no roll design and easy grip handle.
ZEP-46.007.00



Drop Control Blade Holder

Push button blade removal, faster and safer.
ZEP-46.007.02

Miniature Blades

Precise stainless steel mini blades. Individually packed.
Boxes of 25. Sterile.

- SAL-007* Single edge, No.64 round end
- SAL-008* Single edge, No.65 pointed end
- SAL-010* Single edge, No.67 pointed end
- SAL-011* Double edge, No.63 pointed end
- SAL-012* Double edge, No.69 round end



Standard Blades

Stainless steel blades. Individually packed. Boxes of 100. Sterile.



Parallel Double Blade Holder

For cutting tissue strips.
ZEP-46.007.40 (1.5mm strip)
ZEP-46.007.45 (2.0mm strip)



Pivoting Blade Holder

ZEP-46.007.05 Handle
ZEP-46.007.10 Handle and pivot head, right angled



Standard Pivot Head

For above holder.
ZEP-46.007.20



Pivot Head - right angled

ZEP-46.007.50

Bone Scrapers

Safescraper Bone Collection System


At last, a cost effective single use bone scraper!
Single use system for quick and efficient
collection of cortical bone shavings.
Excellent method for collecting autogenous bone
and increasing graft volume.

MX-0702 Straight (Single)

MX-0702-3 Straight (3 Pack)

MX-3987 Curved (Single)

MX-3987-3 Curved (3 Pack)




Micross Mini Bone Scraper

Minimally invasive device for cortical bone removal.
For access to hard to reach and smaller exposed areas.

MX-4049 Single

MX-4049-3 3 pack



Buser Type Bone Scraper, double ended

Three back action bone scrapers
in different widths, to harvest
autologous bone intra-operatively.

ZEP-24.495.01

3mm & 4mm ends

ZEP-24.495.02


4mm & 5mm ends

ZEP-24.495.03 - most popular

3mm & 5mm ends

ZEP-24.916.20 Arkansas Stone (flat)

ZEP-24.950.00 Sharpening oil



Use the original Arkansas Stone
to keep your blades as sharp as
the day you bought them.

**The cost effective way to collect
autologous bone chips**




Bionik Tunneling Instruments

by Prof. Dr. Hürzeler & Dr Zuhr

Helmut Zepf tunneling instruments are perfect for the preparation of flaps for reconstruction of the alveolar ridge or for root cover with a connective tissue graft.



ZEP-46.040.00
Tunneling Set
1/3 Washtray
Upper Jaw
Lower Jaw
Combination



Fine instruments for minimally invasive surgery.

bionik ZEPF



ZEP-46.040.01 Tunneling Instrument for upper jaw, double ended, two sizes, replaceable tips, 18cm straight.



ZEP-46.040.02 Tunneling Instrument for lower jaw, double ended, two sizes, replaceable tips, 18cm offset.



ZEP-46.035.20 Tunneling Instrument combination, straight and offset replaceable tips.

Tunneling Instruments with a new angle

These instruments have an extra angle to facilitate access in the posterior region and are available in two widths 1.8mm and 2.5mm depending on the extent of the tunnel required.



ZEP-46.040.03 Tunnelling Instrument 1.8mm, 45° angled



ZEP-46.040.04 Tunnelling Instrument 2.5mm, 45° angled


TTS - Titanium Trauma Splint - low profile

<i>TTS-1002</i>	single 52mm x 0.2mm	Gold	<i>TTS-1002-10</i>	pack of 10	Gold
<i>TTS-1004</i>	single 100mm x 0.2mm	Gold	<i>TTS-1004-10</i>	pack of 10	Gold
<i>TTS-1006</i>	single 100mm x 0.2mm	Silver	<i>TTS-1006-10</i>	pack of 10	Silver



Flexible, Strong, Simple

Microsurgery Kit



ZEP-46.007.00	Scalpel handle round
ZEP-46.013.00Z	Miniature Blade Holder
ZEP-26.193.20	Bionik Probe Handle
ZEP-24.455.06	North Carolina Probe 1-15 steps
ZEP-46.035.05	Papilex
ZEP-41.864.40	Buser Periosteal Elevator 4.2/3.8
ZEP-41.878.11	Prichard Elevator
ZEP-41.010.17TC	Castroviejo Needle Holder - straight
ZEP-22.814.17	Micro Suture Forceps
ZEP-22.815.17	Cooley Tissue Forceps - straight, serrated
ZEP-46.081.16SC	Joseph Scissors
ZEP-85.180.00	Universal Wash Tray 1/1 275x178x24mm

Prosthetic Kit

Helmut Zepf has developed a Prosthetic Kit allowing you to loosen more than 90% of all screws available on the market.

A real problem solver

ZEP-47.830.00 Prosthetic kit



The Zepf Essential Surgery Kit

The Zepf Essential Surgery Kit contains everything for oral surgery, implantology and periodontal procedures.


Contents:
 MEGADUO mouth mirror
 Bone scraper/dissector
 Gingiva knife/Papillex
 Plugger/Spoon
 Prichards elevator/spatula
 Lucas curette 2.5mm
 Universal curette 2.5mm
 Hygienist curette H6/H7
 Standard/Carolina probe
 Micro-needle holder
 Joseph scissors
 Drop Control round scalpel handle
 Micro-Adson, toothed
 Tweezers, ergonomic
 Aspirator, metal



ZEP-24.961.04 (light red magenta) ZEP-24.961.03 (cobalt blue)

GM Implant Starter Set

A comprehensive starter kit for basic implant and oral surgery procedures. Items selected by a panel of experienced clinicians to get the less experienced surgeon off to a flying start. All items available separately.



GMCK-01

Implant surgery kit

Contents:

Storage box plus

ZEP-46.007.00 Scalpel handle

ZEP-41.864.40 Buser elevator

ZEP-50.226.16 De Bakkey forceps

ZEP-46.081.16SC Joseph scissors

ZEP-37.454.00 Maty retractor

ZEP-37.437.17 Cawood-Minnesota retractor modified

ZEP-41.877.23 Seldin elevator

ZEP-22.272.15 Micro-Adson forceps

ZEP-41.246.15TC Needle holders -Swedish pattern (4/0 - 5/0)

ZEP-46.050.11SC Iris supercut scissors

ZEP-003 Periostome set - includes handle and blades P1, P2 and P3



SinusLift Instrument Set Mini For Internal Sinus Lift



Advantages at a glance:

- atraumatic
- no infection risk in comparison to a open OP
- no wound care or healing interruption
- extremely shortened treatment period
- high success rate. Secure in use
- shortened osseointegration

ZEP-41.848.40



Cacaci SinusLift Elevator, Zepf-Line handle, with graduation
ZEP-41.848.43



Benex® SinusLift Elevator, # 1, Zepf-Line handle, with graduation
ZEP-41.848.41



Benex® SinusLift Elevator, # 2, Zepf-Line handle, with graduation
ZEP-41.848.42

Supplementary Instruments



Kirsch Sinus 7 Plugger, double-ended, 17.5 cm, Spoon
8 x 10 mm, Titanium Zepf-Line
ZEP-41.868.07




Depth Gauge, with graduation 8-10-12-14-16-18 mm,
exchangeable, M4 x 0.5
ZEP-31.691.00

Stiller Sinus Lift Set

Dr Stiller developed his set with a view to coping with even the most problematic sinus structures. With a full range of angled heads for most situations these instruments may also be bent when the unexpected is encountered.

ZEP-24.995.01

Complete Set including
1x washbasket with lid
and partition for Sinus Instrument
Set by Dr Stiller.



*Instruments may be gently
bent to improve access.*



All items available separately.




ZEP-19.651.13 (ø1.5mm)

ZEP-19.651.14 (ø3.0mm)
Surgical Aspirator

ZEP-85.251.03

Mixing cup, stainless steel,
with plastic lid, (ø40mm)



ZEP-41.822.01

Stiller Sinus Elevator, blunt, 4.0mm, bendable, double ended.

ZEP-41.822.02

Stiller Sinus Elevator, blunt, 2.5mm, bendable, double ended.

ZEP-41.822.03

Stiller Sinus Elevator, blunt, 3.0mm, bendable, double ended.

ZEP-41.822.04

Stiller Sinus Elevator, universal, 3.0mm, bendable, double ended.

ZEP-41.822.05

Stiller Sinus Elevator, universal, 6.0mm, bendable, double ended.

ZEP-41.822.11

Stiller Sinus Elevator, sharp, 4.0mm, bendable, double ended.

ZEP-41.822.22

Stiller Sinus Elevator, sharp, 2.5mm, bendable, double ended.

ZEP-41.822.33

Stiller Sinus Elevator, sharp, 3.0mm, bendable, double ended.

ZEP-41.868.07


Kirsch Sinus 7 titanium plugger/spoon

Zepf Cassette System

Stainless steel mesh washtrays




Washtray 1/1 with 2 instrument supports and 1 instrument retainer.
Dimensions 275 x 178 x 24mm.
Suitable for larger instruments.
ZEP-85.180.00



Washtray 1/2 with 2 silicone-rubber supports and instrument retainer, 178 x 140 x 24mm
ZEP-85.181.00



Washtray 1/3 with 2 silicone-rubber supports and instrument retainer,
178 x 90 x 24mm
ZEP-85.182.00



Washbasket 1/1 with universal silicone profile for 15 instruments lengthwise
ZEP-85.995.20

Implant-Ex

ZEPF Explantation Pliers

Implant restoration in the jaw is a blessing for patient care. With advancing age, it is unfortunately not always possible to avoid having to remove an implant again.

There are many indications for this. This can be due to a lack of supply of replacement parts or a clinical indication that makes the explantation of the implant necessary. With the new IMPLANT-EX forceps from Zepf, this implant extraction can be carried out easily, both in the lower jaw and in the upper jaw. The mouth of the forceps is designed in such a way that the milled grooves and the selected radius allow safe gripping of the different implant diameters.

This means that in many cases it is certainly not necessary to remove the implants with trephine drills and this implant extraction can be carried out in a much more patient-friendly manner.



Implant-ex
Explantation Pliers
ZEP-47.834.00

Explantation Pliers
lower jaw, stainless, 155mm
ZEP-47.834.01



“Purgo Biomaterials offer perfect results at cost effective prices”



Purgo now offer world leading biomaterials, membranes and sutures for the full range of clinical demands at very competitive prices.

Founded in 1999, Purgo Biologies strives to become one of the leading global companies in oral health care with its focus on safe biomaterials for soft tissue and bone regeneration.

Based on the specialized experiences accumulated by our outstanding research personnel, Purgo Research and Development Center based in Seoul is thriving to become the best in the world, specifically in the expertise of oral biomaterials for soft tissue and bone regeneration. All members in Research and Development Center are pursuing the optimized technical developments with various clinical studies, cooperative research with the governments, clinicians and educational institutions.

The solutions manufactured by Purgo are gaining fame throughout the world and Purgo's solutions are widely accepted by global dentists from more than 30 countries.

Our production site is complying with the most international quality standards and regularly inspected by international agencies. Each production stage of our biologics solutions are controlled from the selection of the raw material to the final product.

Highly purified, safe porcine biomaterial, granules and syringe applications

ORDER 5 OF
ANY CODE AND
GET 1 EXTRA
FREE



THE Graft™ is a natural, porcine bone mineral matrix for bone reconstruction. It is produced by removal of all organic components from porcine bone. Due to its native structure THE Graft™ likens the physical and chemical aspects of mineralized matrix of human bone. When packed into a bone defect, THE Graft™ gradually resorbs and replaced with bone during the healing process. It is available in cancellous granules packaged in vial. THE Graft™ is sterilized using gamma irradiation.

Unique proprietary manufacturing process removes very effectively potential immunogenic organic elements keeping the natural structure of the matrix.

THE Graft™ quality and safety have been scientifically demonstrated with invitro, in-viva studies, large case study reports and international randomized clinical research. Systematic review and meta-analysis are conducted on THE Graft™ worldwide.

THE Graft™ has established its fame throughout the world, both scientifically and clinically, becoming a popular bone regeneration material.

THEGraft Natural Bone Graft

Item No	Size
PUR-BG-A15	The Graft 0.25-1.0mm porcine 0.15g 0.36cc
PUR-BG-A25	The Graft 0.25-1.0mm porcine 0.25g 0.6cc
PUR-BG-A05	The Graft 0.25-1.0mm porcine 0.50g 1.2cc
PUR-BG-A10	The Graft 0.25-1.0mm porcine 1.00g 2.4cc
PUR-BG-A20	The Graft 0.25-1.0mm porcine 2.00g 4.8cc
PUR-BG-B05	The Graft 1.0-2.0mm porcine 0.50g 1.8cc
PUR-BG-B10	The Graft 1.0-2.0mm porcine 1.00g 3.6cc

THEGraft Natural Bone Graft - syringe packs

Item No	Size
PUR-TG-AS25	The Graft 0.25-1.0mm porcine single 0.25cc
PUR-TG-AS05	The Graft 0.25-1.0mm porcine single 0.5cc
PUR-TG-AS10	The Graft 0.25-1.0mm porcine single 1.00cc
PUR-TG-BS25	The Graft 1.0-2.0mm porcine single 0.25cc
PUR-TG-BS05	The Graft 1.0-2.0mm porcine single 0.5cc
PUR-TG-BS10	The Graft 1.0-2.0mm porcine single 1.00cc



Biocompatible and safe natural bone grafting material

ORDER 5 OF
ANY CODE AND
GET 1 EXTRA
FREE



THE Graft™ Collagen block is composed of porcine derived bone mineral matrix from cancellous bone and atelocollagen from porcine tendon. THE Graft™ Collagen is a bone graft intended to fill, augment, and/or reconstruct periodontal, oral, and maxillofacial defects.

THE Graft™ Collagen bone mineral matrix is similar to physical and chemical aspects of human bone mineralized matrix. Hydrated collagen components have viscosity that facilitate for blending bone mineral matrix. With this characterization, THE Graft™ Collagen can be trimmed and/or molded to the various defect shapes and can be fixed in bone defect site.

As time passes, THE Graft™ Collagen is partially transformed by the osteoclast and osteoblast cells.



Characteristics of THE Graft™ Collagen



1 Easily Moldable

Made of THE Graft™ granules and collagen,

THE Graft™ Collagen is easier to mold than THE Graft™ granules alone. Therefore, THE Graft™ Collagen has better handling property compared to that of THE Graft™, making it possible to adapt grafting materials to various shapes of the defect site with more ease.



2 Optimal osteoconductivity

While retaining better handling properties, THE Graft™ Collagen is able to form sufficient osseous tissue for implant placement and maintain natural volume and great adhesion property which lead to minimum chair time.



3 Predictable clinical results

With great hydrophilicity, THE Graft™ Collagen stabilizes the clot and aids in revascularization of the grafting material in the defect area to increase cell migration efficiency to the mineral substrate. As a result, fast bone formation can be expected, as well as a predictable clinical result.



It is advantageous for shape and space maintenance

THE Graft™ Collagen Block Type

Item No	Size / Volume	
PUR-TCB-01	7x7x7mm	0.34cc
PUR-TCB-02	8x9x10mm	0.72cc
PUR-TCB-03	10x11x12mm	1.32cc

THE Graft™ Collagen Ring Type

Item No	Size / Volume	
PUR-TCR-01	5x3x8mm	0.22cc
PUR-TCR-02	5x4x9mm	0.26cc
PUR-TCR-03	5x5x10mm	0.29cc

Resorbable collagen membrane with high biocompatibility

THE Cover™




ORDER 5 OF
ANY CODE AND
GET 1 EXTRA
FREE


THE Cover™ Flex

THE Cover™ Stiff


THE Cover™ is a resorbable collagen membrane consisting of porcine-originated Type I Collagen.



30x40mm



25x30mm



15x20mm

✓ THE Cover™ is a resorbable membrane originated from a pure Type I Collagen, and it's highly biocompatible and cell-friendly since chemical crosslinking agent was not used.

✓ There are two types of THE Cover™ ;
THE Cover™ Flex offering easy handling experiences and THE Cover™ Stiff more favorable for space maintenance.

THE Cover™ Flex	
Item No	Unit Size
PUR-BP3-1520	15mm x 20mm
PUR-BP3-2530	25mm x 30mm
PUR-BP3-3040	30mm x 40mm

THE Cover™ Stiff	
Item No	Unit Size
PUR-BP5-1520	15mm x 20mm
PUR-BP5-2530	25mm x 30mm
PUR-BP5-3040	30mm x 40mm



Superior handling PTFE non resorbable sutures

Biotex® Non-Resorbable PTFE Suture is comprised of a single-arm, non-resorbable monofilament suture with a stainless-steel surgical needle connected to the suture. The suture is uncoated, undyed and sterile for single use only, composed of 100% PTFE.

- ✓ SOFT HANDLING
- ✓ BIOLOGICALLY INERT
- ✓ NO TANGLE
- ✓ EASY KNOTTING



ORDER 5 OF ANY CODE AND GET 1 EXTRA FREE

Indications

- Bone grafting procedures
- Periodontal surgery
- Guided tissue regeneration
- Ridge augmentation
- Implant surgery
- Soft tissue grafts

Biotex

Item No	USP Size	Length (cm)	Needle Length (mm)	Circle	Point Type
PUR-BT3016	3-0	45	16	3/8	▽
PUR-BT3019	3-0	45	19	3/8	▽
PUR-BT4013	4-0	45	13	3/8	▽
PUR-BT4016	4-0	45	16	3/8	▽
PUR-BT4019	4-0	45	19	3/8	▽
PUR-BT5013	5-0	45	13	3/8	▽
PUR-BT5016	5-0	45	16	3/8	▽
PUR-BTP4013	4-0	45	13	1/2	⊙

Benefits

- Soft and comfortable for patients
- Soft texture for patient comfort
- Reliable closure period
- Superior handling: provides flexibility in the positioning of a square knot. Easy to tie - Easy to remove
- Nonwicking: Elimination of bacterial wicking usually associated to monofilament
- Maintains tensile strength
- PFOA free

Cost effective PTFE membrane for all your GBR procedures

OpenTex® Non-Resorbable PTFE Membrane is a pure medical-grade polytetrafluoroethylene (PTFE) sheet with inert biological features and predictable barrier effect. Due to the smooth surface and small pore size, OpenTex® PTFE Membrane resists the incorporation of bacteria into its structure and eases the removal of the membrane.

Non-resorbable membrane is sustainable for surgical procedure with no primary closure. OpenTex® Membrane is ideal for space-making feature providing enough space for host cells to adhere to grafting materials. OpenTex® is supplied sterile for single use only and available in various sizes.

OpenTex®



ORDER 5 OF ANY CODE AND GET 1 EXTRA FREE

Indications

GBR (Guided Bone Regeneration)

- Simultaneous use of GBR membrane and implants.
- Augmentation around implant placed in immediate extraction sites or delayed extraction sockets.


GBR (Guided Tissue Regeneration)

- Filling of bone defects after root resection, removal of cysts, and removal of retained teeth.




OpenTex	
Item No	Size
PUR-OpenTex_01	24mm x 30mm
PUR-OpenTex_02	17mm x 25mm

Primary Closure



Non-Primary Closure



OpenTex® Main Features



Non-Resorbable

- 100% medical grade PTFE membrane.
- Biologically inert and chemically non-reactive.
- Healing procedure is not interfered with membrane absorption.



Microporous

- Promote the gingival tissue attachment.
- Enhances ease in the interstitial fluid circulation.
- Resist the bacteria infection and fibroblast cells migration.



Minimally Invasive

- Rapid recovery of soft tissue.
- Primary Closure is not necessary.
- Virtually impervious to bacteria.
- Minimum flap reflection or dissection. Safe from bacteria infection, even in the event of the exposure.



Withstands Exposure

- Protect the tissue regeneration site.
- Regenerated underlying tissue can be evaluated.
- Provide a proper environment for the growth of blood vessel and osteogenic cells.

OpenTex® Benefits



Soft Tissue Obtaining



Aesthetic Implant Restoration



Natural Saliva Passage



Minimally Invasive

OpenTex® Strenghts

1 Stability

Non-resorbable PTFE Membrane offers enough healing time to bone regenerative process.

2 Biologically inert

PTFE is soft tissue friendly so it is ideal material as a barrier for bone regenerative process.

3 Withstands to exposure

PTFE Membrane withstands to exposure since it is impervious to bacteria due to their barrier function.

Characteristics of OpenTex®

Most of Oral Bacteria is larger than 1 μm . Open Tex® is micro-porous material that has the pore size small enough to prevent bacterial infiltration. Biocompatible, OpenTex® facilitates cell adhesion on the surfaces. Test performed shows that the surface of Open Tex® is not toxic causing cells to adhere well on the surface.

General Oral Bacteria



> 1 μm

The matter is **PORE SIZE**

Cost effective Titanium reinforced PTFE membranes

OpenTex®-TR

Membrane is composed of 100% polytetrafluoroethylene (PTFE) sheet and grade 1 titanium frame, which are biologically inert and tissue compatible.

OpenTex®-TR Non Resorbable PTFE Membrane with titanium frame is designed to have a suitable surface structure and porosity to prevent integration and passage of bacteria within the interstices of the material, while maintaining space for host cells adhesion to the device.

OpenTex®-TR provides a favorable environment for neovascularization and healing of defects, through repopulating the bone derived cells and protecting the bony defects from migration of the gingival tissue derived cells.

Since the adequate space maintenance is critical to this procedure, the membrane is sufficiently stiff to prevent spontaneous collapse, but also flexible enough to easily conform to tissue contours and reduce perforations of overlying soft tissue.



**ORDER 5 OF
ANY CODE AND
GET 1 EXTRA
FREE**


Indications

- Extraction socket reconstruction
- Bone regeneration
- Where primary closure isn't possible

OpenTex-TR

Item No	Size
PUR-OpenTex-TR_P01	27mm x 25mm
PUR-OpenTex-TR_P02	24mm x 30mm
PUR-OpenTex-TR_P03	17mm x 25mm
PUR-OpenTex-TR_P05	12mm x 24mm
PUR-OpenTex-TR_P06	14mm x 24mm
PUR-OpenTex-TR_P07	30mm x 40mm

Primary Closure



OpenTex® Main Features

OpenTex®-TR



Non-Resorbable



Minimally Invasive



**Optimal rigidity for
space maintenance**



17x25mm
Open Tex-TR_P01



24x30mm
Open Tex-TR_P02



17x25mm
Open Tex-TR_P03



12x24mm
Open Tex-TR_P05



14x24mm
Open Tex-TR_P06



30x40mm
Open Tex-TR_P07

OpenTex®-TR Benefits

1 Optimal rigidity and strength for space making

OpenTex®-TR is optimal product which is able to be trimmed easily and it is solid enough for space making since it is reinforced with titanium frame.

2 Diverse embedded titanium frame

OpenTex®-TR is designed in various shapes to meet surgeon's demand.

3 Excellent tissue interaction

Its micro porous structure helps the tissue interaction.

3 Easy of use


OpenTex®-TR can be trimmed easily and also removed easily.

- ✓ Membrane can be molded and shaped for tenting and space maintenance.
- ✓ The rigidity of the membrane is enhanced to be used for space maintenance.
- ✓ Provides additional stability in large, non-space-making osseous defects.
- ✓ Provide with little memory of Titanium frame, which enables easy placement of the membrane.
- ✓ Ability to withstand exposure.

**ORDER 5 OF
ANY CODE AND
GET 1 EXTRA
FREE**

Choukroun Quattro PRF System

A new concept in Smart Blood Generation, reduces the need for expensive xenografts and membranes.



RN-APRF 3020 Simple System

Includes centrifuge, steribox, PRF preparation box, instrument set, 20 blood collectors, 100 A-PRF (red) tubes, 50 i-PRF tubes (orange).

RN-APRF 3010 Full System

As above with additional steribox, PRF preparation box and instrument set.

The Concept

Centrifugation of whole blood without anticoagulants to produce materials which accelerate bone and soft tissue healing by

- the action of fibrin as a provisional extracellular matrix
- stimulating the release of growth factors into the site
- the presence of mesenchymal stem cells

CHOUKROUN PRF COURSE WITH DR. STAVROS ELEFThERIOU

1ST MARCH 2023 LONDON
24TH MAY 2023 MANCHESTER
27TH SEPT 2023 LONDON

COURSE PRICE
£695
INC VAT

PLUS GET


FREE

£600 WORTH OF
PRODUCT, WHEN
ORDERING THE
PRF SYSTEM



VISIT WWW.GENERALMEDICAL.CO.UK TO BOOK YOUR PLACE

The Choukroun Quattro system now offers six preset programmes



A-PRF+ PRF clots for membrane production or for cutting and using as a biomaterial either alone or mixed with xenografts - red tubes.



i-PRF Liquid PRF for injection – green tubes.



i-PRF M Liquid PRF for male blood to compensate for different haematocrit levels – green tubes.



i-PRF+ Super i-PRF in larger volume (13ml) purple tubes for facial aesthetics and orthopaedic injections.



A-PRF liquid Liquid A-PRF for mixing directly with granular materials to make sticky bone – red tubes.



L-PRF The old L-PRF protocol to obtain PRF clots for membranes – red tubes

RN-A-PRF+ Red tubes for A-PRF+, box of 100

RN-i-PRF Green tubes for i-PRF and sticky bone production, box of 24

RN-i-PRF+ Purple tubes (13ml) for i-PRF+ aesthetic and orthopaedic joint injections, box of 24

RN-APF2040 Blood collectors, box of 25



Call or email for full information pack and prices

Make larger membranes or sticky bone - Easily!

RN-FAST

Side support for horizontal and vertical augmentation.



RN-FAST FAST Plating and screw kit, contains Steribox, screwdriver holder, c/a and manual drivers, 2 x Kelly clamps, tweezers, 19 screws from 3.5mm to 11mm, 2 x 5 hole plates, 1 x 12 hole plate, 2 hooks.

RN-MESH5T	Mesh 5 holes
------------------	--------------

RN-MESH12T Mesh 12 holes

RN-MESH23T Mesh 23 holes

RN-Q3.5 Screw 3.5mm (5)

RN-Q4 Screw 4mm (5)

RN-Q5 Screw 5mm (5)

RN-Q6 Screw 6mm (5)

RN-Q7 Screw 7mm (5)

RN-Q8 Screw 8mm (5)

RN-Q9 Screw 9mm (5)

RN-Q11 Screw 11mm (5)



Screw 2mm in lengths 3.5, 4, 5, 6, 7, 8, 9, 11mm

All parts available separately

Soft Brushing kit

RN-SBK

Gently extends the periosteum for easier graft site coverage without cutting.



Wen Tension Release Comb

ZEP-41.865.01

Width of working tips 5 and 8mm



Monoglyc suture range

- RN-JCG65116** Monoglyc resorbable monofilament, 4/0 DS18 70cm Violet box 24
- RN-JCG61417** Monoglyc resorbable monofilament, 4/0 DS21 70cm Violet box 24
- RN-JCG65120** Monoglyc resorbable monofilament, 4/0 DS18 45cm Violet box 24
- RN-JCG65121** Monoglyc resorbable monofilament, 5/0 DS18 45cm Violet box 24
- RN-JCG61507** Monoglyc resorbable monofilament, 5/0 DSM13 70cm Violet box 24
- RN-JCG65109** Monoglyc resorbable monofilament, 5/0 DSM16 45cm Violet box 24
- RN-JCG65122** Monoglyc resorbable monofilament, 5/0 DS18 70cm Violet box 24
- RN-JCG65127** Monoglyc resorbable monofilament, 6/0 DSM13 45cm Violet box 24

Osteo-Mesh™ TM-300 Titanium Mesh

- Non-stick nitride coating
- Unused portions can be re-sterilized by autoclave
- Porosity of 0.5 mm contains most graft particles
- Strong titanium construction allows for space-maintenance

CYT-OSTEOMESHTM25X34

Osteo-Mesh TM-300 25mm x 34mm

CYT-OSTEOMESHTM45X45

Osteo-Mesh TM-300 45mm x 45mm

**Best secured by using Master Pin - Page 12
or Pro-fix Kit - Page 18**



Tidespon

Absorbable haemostatic gelatin sponge, made up of 100% porcine gelatin

Alternative to our best selling collagen cubes

Resorbable haemostatic sponge cube. Manufactured from 100% porcine gelatine. Adheres to the wound/extraction site absorbing 50% it's weight. Activates the clotting cascade and prevents further bleeding while providing a focus for regeneration, completely resorbs in a few weeks.

AS-TS-10-10-10 Tidespon Dental Sponge, 10x10x10mm, box of 50

NEW



Cytoplast RTM Wound Dressings

The Cytoplast® RTMPlug, RTM Tape, and RTMFoam absorbable wound dressings consist of highly purified collagen that is derived from bovine tissue of United States origin. The RTM wound dressings are essentially resorbed within 30 days. The wound dressings are used in the treatment of surgical wounds, extraction sites, dental sores (non-infected or viral) and oral ulcers.

CYT-RTMPLUG10 Cytoplast RTMPlug, 1 x 2cm, pack of 10

CYT-RTMTAPE10 Cytoplast RTMTape, 2 x 4cm, pack of 10

CYT-RTMFOAM10 Cytoplast RTMFoam, 2.5 x 7.5cm, pack of 10

NEW



Cytoplast GBR-200




The original, patented, high density PTFE membrane.

- Medical grade dense PTFE membrane
- Non-resorbable, yet easily removed without surgery
- Primary closure not required <0.2um pore size; virtually impervious to bacteria
- No risk of bacterial invasion due to exposure of membrane
- No disruption of healing due to resorption process
- Biologically inert and chemically non-reactive 12 x 24mm size is ideal for extraction sites

<i>CYT-GBR-200-S</i>	GBR-200 membrane, single 12 x 24mm
<i>CYT-GBR-200-S-10</i>	GBR-200 membrane, box of 10 singles 12 x 24mm
<i>CYT-GBR-200-L</i>	GBR-200 membrane, single 25 x 30mm
<i>CYT-GBR-200-L-4</i>	GBR-200 membrane, box of 4 singles 25 x 30mm

**ADD OUR
TACKSYSTEM**





Cytoplast TXT-200

Textured dense PTFE membranes.

The textured surface increases surface area for greater membrane stabilisation.
Withstands exposure and makes socket preservation more predictable.

<i>CYT-TXT1224</i>	Textured membrane, 12 x 24mm, pack of 10
<i>CYT-TXT1224-1</i>	Textured membrane, 12 x 24mm, single
<i>CYT-TXT2530</i>	Textured membrane, 25 x 30mm, pack of 4
<i>CYT-TXT2530-1</i>	Textured membrane, 25 x 30mm, single





RPM Reinforced PTFE mesh

RPM's unique macroporous design allows for direct contact between the bone graft and periosteum, allowing naturally occurring revascularization and infiltration of cells into the bone graft.

Versatile Rectangular Shapes



PSE

CYT-RPM200PSE
20mm x 25mm



PLE

CYT-RPM200PLE
25mm x 30mm



XLE

CYT-RPM200XLE
30mm x 40mm



XLKE

CYT-RPM200XLKE
30mm x 40mm



XLKME

CYT-RPM200XLKME
30mm x 40mm



K2E

CYT-RPM200K2E
40mm x 50mm

Shapes with Fixation Points



BLE

CYT-RPM200BLE
17mm x 25mm



PSTE

CYT-RPM200PSTE
25mm x 36mm



PLTE

CYT-RPM200PLTE
30mm x 41mm

Interproximal Shapes



ATCE

CYT-RPM200ATCE
24mm x 38mm



ATCME

CYT-RPM200ATCME
24mm x 38mm



PTCE

CYT-RPM200PTCE
38mm x 38mm



PTCM

CYT-RPM200PTCM
38mm x 38mm



PDE

CYT-RPM200PDE
38mm x 38mm



PDMR

CYT-RPM200PDMR
38mm x 38mm



PDML

CYT-RPM200PDML
38mm x 38mm

CytoPlast Dense PTFE Titanium Reinforced Membranes



New 150mm micron thickness

Rectangular



CYT-Ti150ANL
12mm x 24mm



CYT-Ti150PS
20mm x 25mm



CYT-Ti150PL
25mm x 30mm



CYT-Ti150XL
30mm x 40mm



CYT-Ti150XLK
30mm x 40mm



CYT-Ti150K2
40mm x 50mm



CYT-Ti150AS
14mm x 24mm

Original 250mm micron thickness



CYT-Ti250ANL
12mm x 24mm



CYT-Ti250ANL30
12mm x 30mm



CYT-Ti250PS
20mm x 25mm



CYT-Ti250PL
25mm x 30mm



CYT-Ti250XL
30mm x 40mm



CYT-Ti250XLK
30mm x 40mm



CYT-Ti250K2
40mm x 50mm



CYT-Ti250AS
14mm x 24mm



CYT-Ti250ATC
24mm x 38mm

**BEST
SECURED
USING PRO-FIX**
SEE PAGE 18

Cytoplast PTFE - the soft monofilament suture

Cytoplast™ PTFE suture is an ideal suture for dental bone grafting and implant procedures where a soft monofilament suture is desirable. Cytoplast™ PTFE suture's monofilament construction doesn't allow bacterial wicking into the surgical site, and patients like it because it is soft.



100% medical grade PTFE

Biologically inert, chemically non reactive, superior soft tissue response.

Soft monofilament

No bacterial "wicking" and comfortable for patients.

Non absorbable

Reliable closure period.

No package memory

Excellent handling and secure knots.

<i>CYT-CS-04</i>	18" 2/0 19mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-0428</i>	28" 2/0 19mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-05</i>	18" 3/0 16mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-06</i>	18" 4/0 13mm 3/8 premium reverse cut. (pack of 12)
<i>CYT-CS-0628</i>	28" 4/0 13mm 3/8 premium reverse cut. (pack of 12)
<i>CYT-CS-07</i>	18" 3/0 19mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-07 Black*</i>	18" 3/0 19mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-0728 Black*</i>	28" 3/0 19mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-08</i>	18" 4/0 16mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-0828</i>	28" 4/0 16mm 3/8 reverse cut (pack of 12) (for delicate tissues requiring atraumatic needle)
<i>CYT-CS-10</i>	18" 5/0 13mm 3/8 reverse cut (pack of 12)
<i>CYT-CS-11</i>	18" 5/0 16mm 3/8 reverse cut (pack of 12)

* Black needles provide more contrast in the mouth for easier visualisation

Sutures

Most
Popular
Sizes

PROLENE (Substitute for Seralon)

Monofilament, non resorbable

CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
W8695H	13mm	6/0	45cm	24
W8868T	11mm	6/0	45cm	24
W8882T	16mm	5/0	45cm	24
W8872T	13mm	5/0	45cm	24
W8683	19mm	4/0	45cm	12
W8884T	16mm	4/0	45cm	24
W8005T	16mm	6/0	45cm	24

VICRYL *rapide*

Braided resorbable half strength 5 days

CODE	NEEDLE	SIZE	LENGTH m	PACK SIZE
W9974	3/0	3/0	0.45	15
W9913	11 3/8 P	6/0	0.45	12
W9915	11 3/8 P	5/0	0.45	12
W9910	16 Compound	5/0	0.45	12
W9911	16 Compound	4/0	0.45	12
W9918	16 3/8 P	4/0	0.75	12
W9922	19 3/8 P	4/0	0.75	12
W9930	19 3/8	4/0	0.75	12
W9935	22 1/2	3/0	0.45	12
W9951	19 3/0	4/0	0.45	36

MONOCRYL

Monofilament resorbable half strength 7 days

CODE	NEEDLE	SIZE	LENGTH m	PACK SIZE
W3214	11 3/8 P	6/0	0.45	12
W3215	13 3/8 P	6/0	0.45	12
W3209	19 3/8	5/0	0.70	12
W3206	19 3/8 P	4/0	0.45	12
W3207	19 3/8 P	3/0	0.45	12

Wide range of smaller and larger size sutures available on request

Coated VICRYL

Braided resorbable half strength 21 days

CODE	NEEDLE	SIZE	LENGTH m	PACK SIZE
W9500T	11 3/8 P	6/0	0.45	24
W9501T	11 3/8 P	5/0	0.45	24
W9436	16 compound	5/0	0.45	12
W9506T	16 3/8 P	4/0	0.45	24
W9437	16 compound	4/0	0.45	12
W9386	19 3/8	4/0 violet	0.45	12
W9951	19 3/8	4/0	0.45	12
W9471	22 3/8	4/0 violet	0.45	12
W9472	22 1/2	3/0 violet	0.45	12
W9741	22 1/2	3/0	0.45	12
W9824	17 3/8	3/0	0.75	12

ASSUPRO (SERALON/PROLENE type)

Non resorbable polypropylene monofilament

		NEW		BEST VALUE
CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
AS-PY-813	16mm 3/8 Precision needle	6/0	45cm	36
AS-PY488	13mm 3/8 Precision needle	6/0	45cm	36
AS-PX489	13mm 3/8 Precision needle	5/0	45cm	36
AS-PX817	16mm 3/8 Precision needle	5/0	45cm	36
AS-PY485	13mm 3/8 Standard needle	6/0	45cm	36
AS-PY498	16mm 3/8 Standard needle	6/0	45cm	36
AS-PX503	13mm 3/8 Standard needle	5/0	45cm	36
AS-PX499	16mm 3/8 Standard needle	5/0	45cm	36
AS-PW515	18mm 3/8 Standard needle	4/0	45cm	36
AS-PW810	13mm 3/8 Precision needle	4/0	45cm	36
AS-PW484	11mm 3/8 Precision needle	6/0	45cm	36

Wide range of smaller and larger size sutures available on request

ASSUNYL POLYAMIDE

Non resorbable polyamide monofilament

CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
YY485	13mm 3/8 reverse standard needle	6/0	45cm	36
YX503	13mm 3/8 reverse standard needle	5/0	45cm	36
YX499	16mm 3/8 reverse standard needle	5/0	45cm	36
YW500	16mm 3/8 reverse standard needle	4/0	45cm	36
AS-YY813	11mm 3/8 reverse precision needle	5/0	45cm	36

GLICOFIL LAC (SERAFIT/VICRYL COATED type)

Multifilament resorbable, half strength 15-20 days

CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
AS-FW870LW	16mm 3/8 Precision needle	4/0	45cm	36
AS-FX503L	13mm 3/8 Standard needle	5/0	45cm	36
AS-FX499L	16mm 3/8 Standard needle	5/0	45cm	36
AS-FW500L	16mm 3/8 Standard needle	4/0	45cm	36
AS-FX489LW	17mm 3/8 Precision needle	5/0	45cm	36

GLICOFIL LAC fast (SERAPID/VICRYL RAPIDE type)

Fast resorbing multifilament, half strength 5-7 days

CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
AS-FX511LRAH	16mm 3/8 Standard needle	5/0	45cm	36
AS-FW812LR	18mm 3/8 Standard needle	4/0	45cm	36
AS-FY484LR	11mm 3/8 Precision needle	6/0	45cm	36

MONOFIL (SERAFAST/MONOCRYL type)

Resorbable monofilament

CODE	NEEDLE	SIZE	LENGTH cm	PACK SIZE
AS-FX489FW	13mm 3/8 Precision needle	5/0	45cm	36
AS-FX517FW	18mm 3/8 Standard needle	5/0	70cm	36
AS-FW500FW	16mm 3/8 Standard needle	4/0	70cm	36

Wide range of smaller and larger size sutures available on request

**NEW**

Narrow Ridges? Problem Solved!

Why should you be forced to use round implants when other shapes would better fit your patients' anatomy?





A new shape in implantology with
REX PiezoImplant




The obvious choice?

mectron

Advantages of Wedge Implants



Implant site preparation performed with piezoelectric tips



Wedge implants seating in the final position.



Globally Trusted, Now for FIXED

The Gold Standard for Standard and Narrow Ridges

The LOCATOR® Connection

The LOCATOR® implant system, now supporting both removable and FIXED treatment plans, offers you the confidence that comes with one of the most trusted names in dentistry. LOCATOR and LOCATOR FIXED use the same LOCATOR abutment and workflow, which allows you to expand your patient universe and grow your practice. By bundling the abutment, processing components, and implant we're able to deliver an implant system that will change the way you approach edentulous or partially edentulous patients. This all-in pricing can reduce the cost of LOCATOR implant and abutment pairings by 50% or more, increasing the universe of treatable patients.

Best of all, they come in the diameters needed for both standard and narrow ridged patients.

Available Diameter sizes: 2.4, 2.9, 3.5, 3.9, 4.4, 4.9mm

LOCATOR Removable



LOCATOR FIXED



2.4mm

In addition to simplifying ordering and inventory, LOCATOR's all-in approach allows you to provide full-arch solutions to a range of underserved patients. The patient-friendly pricing can allow you the option of doing 'four implant overdenture procedures' at the same cost of doing a 'two-implant procedure' with other higher priced implants or reduce the cost of removable overdenture and Fixed full-arch procedures to reach more potential patients. All while maintaining the biomechanical standards that have earned the trust of dental professionals.



FOR NARROW & STANDARD RIDGE CASES



The LOCATOR® Advantages

**ZD ZEST DENTAL
SOLUTIONS®**

RBM

Surface treatment for high level of bone to implant contact and osseointegration

Tapered Core

Gradual bone condensing for reliable functionality

Tip Design

Custom tip design based on implant diameter for ease of placement

One Connection for All

Single 2.9mm prosthetic connection for all implant sizes enabling a platform switch for standard ridge designs

Dual Lead Custom Thread

Progressive thread design for primary bone stability, bone condensing and speed of insertion. Thread design is more aggressive as implant diameter increases.



Driven by Innovation



NEW LOCATOR Implant System is compatible with LOCATOR R-Tx Abutments offering 25%1 greater wear resistance, an impressive 60° of total implant divergence and dual retentive surfaces for easier seating

of the overdenture. LOCATOR R-Tx delivers greater durability with easier maintenance and is now the perfect companion to the AFFORDABLE LOCATOR Overdenture Implant!

Easy-to-Follow Surgical Kits



From Removable to **FIXED**... Upgrade to LOCATOR **FIXED** with the Same Implant and Abutment

Introducing a NEW flexible and accessible **FIXED** full-arch solution

- LOCATOR Removable or **FIXED**... it's your choice.
- Utilize patients' existing or newly placed LOCATOR Abutments to transition to a **FIXED** full-arch solution.
- Offer staged treatment planning for your patients, removable today... **FIXED** tomorrow.
- Give your patients the life they had with their permanent teeth with an affordable, permanent teeth solution
- A LOCATOR proven alternative to the discomfort and pain associated with conventional dentures that feels more like permanent teeth.



**Same Abutment,
Reduced Chairtime,
Less Bone Reduction...
LOCATOR **FIXED****

Legacy Locator Implant Abutment System

The original Zest Locator system used for more than 15 years by leading implantologists worldwide. Fits all major implant systems.

LOCATOR ABUTMENTS/LASER BAR

IMPLANT ABUTMENTS

Implant Abutment
Multiple Implant Systems and
Tissue Cuff Heights



Multi-Unit Abutment
w/ Titanium Collar



Bar Abutment Thread
08587 2-56
08589 2.0mm



Cast-to Abutment
Stainless Steel
08586



ROOT ABUTMENTS

0° Root Abutment
Stainless Steel
08520



10° Root Abutment
Stainless Steel
08521



20° Root Abutment
Stainless Steel
08522



Cast-to Coping
Stainless Steel
08528



LASER BAR

Female Stainless Steel
08590-02
08590-10



Female Titanium
08588-02
08588-10



LAB PROCESSING

Impression Coping
08505



Abutment Analog
08530 4mm
08516 5mm



Abutment
w/Delrin Collar*
.08917



Block-Out Spacer
08514



White

Processing Spacer
08569



Parallel Post
08517



Denture Cap Assembly
08510



Standard Range
Male Processing Pkg
08519



Extended Range
Male Processing Pkg
08540



Bar Male
Processing Pkg
08028



Black
Processing Male
08515



Black

Yellow Bar
Processing Male
08026



Yellow

Castable Threaded Insert
08013 2-56
08014 2.0mm



Delrin

REPLACEMENT MALES

STANDARD RANGE

Extra Light
Retention Male
08529



Blue

Light
Retention Male
08527



Pink

Regular
Retention Male
08524



Clear

Zero
Retention Male
08558



Grey

Extra Light
Retention Male
08548



Red

Light
Retention Male
08915



Orange

Regular
Retention Male
08547



Green

DRILLS & TAPS

Spot Face
Diamond Bur
08922



Pilot Drill
08924



Bar Drill
09102 1.7mm
(2.0mm Thread)
09103 1.8mm
(2-56 Thread)



Bar Tap
09104 2.0mm
09105 2-56



Drill & Tap
Holder
08016



Paralleling Mandrel
09107



Angle
Measurement Guide
09530



TORQUE DRIVERS & TOOLS

Torque Wrench Kit:
Torque Wrench, 15mm Square
Drive Insert and Thumb Knob
04391 20Ncm
9020 30Ncm



Square Drive Torque
Wrench Driver
08926 15mm
08927 21mm



Latch Type Torque
Wrench Driver
08913 23mm
08914 29mm



LOCATOR Core Tool:
Male Removal Tool, Male Seating Tool
& TiN Coated Abutment Driver
08393



LOCATOR **FIXED**™

ZD ZEST DENTAL SOLUTIONS®



Give your patients the life they had with their permanent teeth with a less invasive and less expensive solution.


- No New Abutment Required...
- Patient can easily transition between fixed and removable solutions without changing abutment
- Increase case acceptance for patients who do not want a denture but cannot afford traditional screw retained solutions
- Chairtime savings due to no screws or cement needed for installation or hygiene appointments (prosthesis can be removed and seated in seconds)
- Less need for bone reduction due to shorter prosthetic space requirements vs other fixed systems

UPGRADE TO LOCATOR **FIXED**

LOCATOR FIXED Prosthesis Housing is designed to have higher bond strength and maintains the same pick-up method as existing LOCATOR housings.

Ultra precision manufactured FIXED inserts made of rigid material designed to match the inner geometry of our LOCATOR FIXED housing performs as a fixed full arch solution. Corrects up to 40 degrees of total implant divergence.

LOCATOR Abutment - Build on what you already know and trust in your practice. Zest LOCATOR FIXED™ delivers a fixed full arch solution using existing LOCATOR abutments.



Maintains same prosthetic space requirements as LOCATOR Removable System!

New LOCATOR **FIXED**...Same Workflow as LOCATOR Removable

Enhanced LOCATOR System Core Tool

One streamlined tool compatible with LOCATOR removable and LOCATOR FIXED Inserts.

Two-sided instrument designed for easy insertion and removal of any LOCATOR Insert.



Locator R-Tx


Better. Simpler. Stronger.

ZD ZEST DENTAL
SOLUTIONS®



Unveiling LOCATOR R-Tx, a better, simpler, and stronger system that relies on the same restorative techniques as the award-winning LOCATOR.

- Based on 40 years experience with the undisputed market leader for implant retained overdentures.
- Industry standard .050 hex drive placement.
- Less food and plaque accumulating.
- Stronger and harder wearing for longer life.
- Aesthetic gingiva coloured abutment.



CUSTOM RETENTION INSERT TOOL

















Dual-sided tool for easy insertion and removal of the LOCATOR R-Tx Retention Insert. Metal slot feature on the back of the Tool allows the clinician to disengage a Retention Insert without the risk of tearing gloves.



CONVENIENT ALL-IN-ONE PACKAGING

Custom designed all-in-one, double-ended vial separately holds Abutment and Processing Components providing all the necessary components for the case with one part number.

Locator R-Tx

CODE No.	DESCRIPTION	
PROCESSING COMPONENTS		
LOC-30001	Zero Retention Insert (Qty. 4)	
LOC-30002	Low Retention Insert (Qty. 4)	
LOC-30003	Medium Retention Insert (Qty. 4)	
LOC-30004	High Retention Insert (Qty. 4)	
LOC-30012	Processing Insert (Qty. 4)	
LOC-08514	Block-Out Spacer (Qty. 20)	
LOC-30013	Denture Attachment Processing Assembly (Qty. 4)	
LOC-30014	3.35mm Abutment Analog (Qty. 4)	
LOC-30015	4mm Abutment Analog (Qty. 4)	
LOC-30016	5mm Abutment Analog (Qty. 4)	
LOC-30017	Impression Coping (Qty. 4)	
LOC-30018	Processing Spacer (Qty. 4)	
TOOLS		
LOC-30021	Locator R-Tx Retention Insert Tool (Dual Sided)	
LOC-08007	.050 "/1.25mm Hex Driver with thumb knob (Shaft 12.7mm)	
LOC-08008	.050 "/1.25mm Hex, Latchlock Driver (Shaft 12.7mm)	
LOC-04936	.050 "/1.25mm Hex, Torque Wrench Insert (Shaft 6.4mm)	

Chairside APM



The CHAIRSIDE attachment processing material you've been waiting for.

A tissue-colored, dual-cure and self-adhesive composite indicated for the processing of attachments into full and partial overdentures, using either a chairside or a laboratory procedure. Available in a 8gm syringe and 18mL cartridge, the Attachment Processing Material enables easy pick-up of all LOCATOR Denture Attachment Housings/Caps.

Simple

- No primer needed for reduced chairtime
- Bonds to itself for easy void-fills and touch-ups
- Dual-cure capability provides convenience and flexibility

Versatile

- Normal Set and Fast Set varieties match the complexity and speed of your procedure
- Available in 8gm syringes or 18mL cartridges
- Shorter mixing tips and ideal viscosity for precise application, as well as reduced

waste and cost

Patient Friendly

- Low curing temperature
- Odorless and tasteless

LOC-9566 Chairside 4ml (8gm) syringe and 15 mixing and angled tips

LOC-9565 Chairside 18ml syringe with 10 blue mixing tips

LOC-9125-15 Short mixing tips for use with 8gm syringe



Chairside Fast set

A work time of 30 seconds and a set time of 2 minutes and makes picking up attachments faster and easier.

LOC-90-00130 Chairside Fast Set



Chairside Prep & Polish Kit

- A comprehensive kit consisting of Recess, Trim, Undercut, Grind and Vent Burs, as well as a Polisher, all of which are designed to address the most frequent overdenture preparation requirements.
- Single use CHAIRSIDE Recess Bur, specifically designed by ZEST Anchors, that quickly and easily prepares the exact size recess for the LOCATOR and SATURNO Denture Attachment Housing (Denture Caps).

LOC-9582 Chairside Denture Prep & Polish Kit

LOC-9580-5 Pack of 5 single use polishers

LOC-9576 Recess Bur

LOC-9577 Undercut Bur

LOC-9578 Vent Bur

LOC-9579 Trim Bur

LOC-9583 Grind Bur



Improv Temporary Implant Cement

- **IMPROV** offers strength & predictability by reducing the risk of loose crowns and improper seating.
- **IMPROV's** dual barrel delivery syringe eliminates inaccuracy and the mess of two tube mixing.
- **IMPROV** offers a consistent working time without premature setting.
- **IMPROV** is Eugenol free and can be effectively removed easily from both hard and soft tissue.
- **IMPROV** is the ideal provisional cement for both implants and natural dentition.



USED
BY THE
EXPERTS

SAL-IMPROV

Contains 2 x 6ml Syringes

12 x Mixing Tips, 6 x spatulas, 80 x Mixing Pads

SAL-IMPLCEMITIPS Extra mixing tips x 30

Locator Fixed Seating and Removal Tool

Locator Fixed Seating and Removal Tool (includes Tool, Seating Tip, Removal Tip with Threaded Adapter, Wire/Lever Wrench and Tip Wrench)

LOC-10100




SaturnO O-Ring Attachment

SaturnO provides the perfect replacement for conventional O-ring fixtures offering greater angle correction (up to 30 degrees) on divergent implants. SaturnO attachments last much longer than other O-rings due to their robust construction and also give a more secure fit to the denture, a feature greatly appreciated by patients.


Offered in Standard (2.25mm ball) and Micro (1.8mm ball) sizes SaturnO are easily fitted to all major O-ring attachments.

SATURNO

**Fits
any ball
system**



- LOC-1966-2** Standard (for 2.25mm ball) attachment assembly, (pack of 2)
- LOC-1966-10** Standard (for 2.25mm ball) attachment assembly, (pack of 10)
- LOC-1969-10** Replacement standard rings (for 2.25mm ball) attachment assembly, (pack of 10)
- LOC-1971-2** Micro (for 1.8mm ball) attachment assembly, (pack of 2)
- LOC-1971-10** Micro (for 1.8mm ball) attachment assembly, (pack of 10)
- LOC-1974-10** Replacement Micro rings (for 1.8mm ball) attachment assembly, (pack of 10)
- LOC-8049** Saturno O-Ring insertion tool



Locator Attachment System for Multi-Unit Abutments

Convert with Confidence


Transition from a screw-retained fixed prosthesis to a removable overdenture solution with ease by utilising the existing implant manufactured multi-unit abutments with the Locator Attachment System for Multi-Unit Abutments.



Self-Retained Solution with existing implant Manufacturer Multi-Unit Abutments



Removable overdenture solution with the Locator Attachment System for Multi-Unit Abutments



Note: 5.35mm of vertical clearance is required from the shoulder of the multi-unit abutment.

- Self-aligning, dual retentive design to enhance engagement and stability of the overdenture
- Low profile to maximise interocclusal space
- Customisable levels of retention

Low cost option for fixed to removable prosthesis

Designed for free-standing and castable bar-splinted applications

Multi-Unit Abutment w/Titanium Collar



Male Processing Pkg



Extended Range Male Processing Pkg



Abutment w/Delrin Collar




Bar Male Processing Pkg



Part Number	Description
Free-Standing Application	
08909-2	LOCATOR Abutment and Collar for Multi-Unit (2 Pack or 10 Pack)
08909-10	
08519-2	LOCATOR Male Processing Package (Includes Denture Cap with Black Processing Male, Standard Range Males: Clear, Pink, Blue and Block-Out Spacer) (2 Pack or 10 Pack)
08519-10	
08540-2	LOCATOR Extended RangeMale Processing Package (Includes Denture Cap with Black Processing Male, Extended Range Males: Green, Orange, Red and Block-Out Spacer) (2 Pack or 10 Pack)
08540-10	
Bar Splinted Application	
08917-2	LOCATOR Abutment and Castable Collar for Multi-Unit (2 Pack or 10 Pack)
08917-10	
08028-2	LOCATOR Bar Male Processing Package (Includes Denture Cap with Yellow Processing Male, Standard Range Males: Clear, Pink, Blue and Block-Out Spacer)
08028-10	

For details of the Locator multi-unit conversion system contact the company.

Sterile Packs for all procedures



GM-BASIC

Sterile procedure pack with all the essentials for implant, oral surgery and advanced perio cases. Quality components at competitive prices, often copied, never bettered.

- | | |
|---|---|
| 1 Wrapping sheet (100 x 100cm) | 2 Drapes (75 x 75cm) with adhesive edge |
| 1 Table cover (75 x 90cm) | 1 Patient U-Drape (160 x 200cm) |
| 2 Face masks | 2 Hose/tubing covers (7 x 120cm) |
| 1 Patient bouffant cap | 2 Light handle covers (15 x 20cm) |
| 1 Surgeons cap with ties | |
| 1 Nurses cap (elasticated) | <i>DH - Basic Set - 1</i> |
| 2 Gowns Classic, large inc 4 towels | <i>DH - Basic Set - 8</i> |
| (Gowns & towels in separate wrap within pack) | <i>DH - Basic Set - 32</i> |

MOST POPULAR

GM-EXPOSURE


Designed for implant exposure and minor surgeries.

- 1 Wrapping sheet (100 X 100cm)
- 2 Prima surgical gowns (LARGE) with towels
- 1 Mini U-drape, 58 X 80cm
- 1 Surface drape with adhesive (75 x 90cm)
- 2 Hose/tubing sleeves
- 3 Swabs (7.5 x 7.5) packs of 5

DH-60-19-170-S Box of 8

Check out our competitive pricing for Basic and Custom Procedure Packs

Sterile Packs for all procedures



GM-COMP

Our own brand comprehensive procedure pack, all you need at an affordable price.

- | | |
|--|---|
| 1 Wrapping sheet/Tyvec pack (100 x 100cm) | 2 Swabs (10 x 10cm) pack of 5 |
| 1 Table cover | 3 Swabs (5 x 5cm) pack of 5 |
| 2 Surgical gown Classic (large) including 4 towels | 1 Yankauer tip small with vac control |
| 1 Patient U-drape (160 x 200cm) | 1 Hygovac/hygosurg combo |
| 2 Surface drape with adhesive (75 x 75cm) | 1 Forceps (plastic) 12cm |
| 2 Light handle cover (15 x 20cm) | 1 Plastic tray - 22 x 15cm - one large 3 smaller sections |
| 4 Hose/Tubing sleeve | 4 Adhesive strips 2.5 x 2.5cm |
| 1 Motor unit cover (70 x 80cm) | 1 Waste bag - yellow |
| 2 Face masks | |
| 1 Bouffant cap | |
| 1 Surgeons cap (ties) | |
| 1 Nurses cap (elasticated at back) | |

GM-COMP-1 Single
GM-COMP-7 Box of 7

GM-UNIVERSAL

For second stage, minor oral and perio procedures.

- 1 Wrapping sheet (100 x 100cm)
- 1 Mini split U-Drape (58 x 80cm)
- 2 Drapes (75 x 90cm) with adhesive edge
- 2 Light handle covers
- 2 Hose/tubing covers

DH -Minor Set - 10

Custom Procedure Packs

Prescription Form

General Medical custom packs are high quality economically priced and made to your personal specification. Each outer pack is made of Tyvec for added security and inside components are presented within a wrapping sheet which acts as an additional table cover. Because of our very low prices we ask that you commit to an initial order of around 16 units. This number is variable depending on the size of each pack and how many can be packed per carton. For a quote request a prescription form or visit our website www.generalmedical.co.uk.

All items on our custom prescription list are available separately and listed in this catalogue.

Product

Wrapping sheet/Tyvec pack (100 x 100cm)
 Table Cover
 Surgical gown classic (medium) + 2 towels
 Surgical gown classic (large) + 2 towels
 Surgical gown classic (extra large) + 2 towels
 Mini U-drape, short procedure
 Patient U-drape (150 x 200cm)
 Patient U-drape (125 x 150cm)
 Patient U-drape (160x200cm) with neck cut out
 Fenestrated drape (10cm window)
 Incise Foil (Tegaderm type) (20 x 30cm)
 Surface drape (75 x 75cm)
 Surface drape with adhesive (75 x 75cm)
 Surface drape (75 x 90cm)
 Surface drape with adhesive (75 x 90cm)
 Surface drape (75 x 120cm)
 Surface drape with adhesive (75 x 120cm)
 Light handle cover (15 x 20cm)
 Plastic Mittens (2)
 Concertina hose cover
 Hose/Tubing sleeve
 Mayo cover (75 x 145cm)
 Motor unit cover (70 x 80cm)
 Motor unit cover elasticated (85 x 85cm)
 Motor cover (15 x 20cm)
 Face masks (ties)
 Face masks (elastic loops)
 Bouffant caps
 Surgeons cap (ties)
 Nurses cap (elasticated at back)
 Swabs (10 x 10cm) pack of 5
 Swabs (5 x 5)cm pack of 5
 Suction Tubing (single) 2m
 Microtip and Tube set 2m
 Yankauer tip small with vac control
 Hygovac/hygosurg Combo
 Forceps (plastic) 12cm
 Plastic tray – kidney bowl type
 Plastic tray – (22 x 15cm) – one large/3 small sections
 Adhesive strip (2.5 x 25cm)
 Gallipots 125ml



**Visit www.generalmedical.co.uk
 for full details of our custom pack
 service**

Gowns

Classic

Three layer hydrophobic non woven gown.
Edged, elasticated sleeves with two towels
included in the pack.

Medium	<i>DH-351023-21-S-25</i> (box of 25)
Medium	<i>DH-351023-21-S</i> (single)
Large	<i>DH-351033-21-S-25</i> (box of 25)
Large	<i>DH-351033-21-S</i> (single)
Extra Large	<i>DH-351043-21-S-25</i> (box of 25)
Extra Large	<i>DH-351043-21-S</i> (single)

Prima

Lightweight spunbond non woven gown for low fluid
procedures, tie tape closure and elasticated cuffs.
Two towels included with each gown.

Medium	<i>DH-354523-20-S-30</i> (box of 30)
Medium	<i>DH-354523-20-S</i> (single)
Large	<i>DH-354533-20-S-30</i> (box of 30)
Large	<i>DH-354533-20-S</i> (single)
Extra Large	<i>DH-354543-20-S-30</i> (box of 30)
Extra Large	<i>DH-354543-20-S</i> (single)



Caps

Surgical Caps

- DH-15520* Male cap with ties, blue, pack of 100
DH-16624 Male cap, elasticated, blue, pack of 100
DH-14522 Female cap with ties, green, pack of 100
DH-17214 Female cap, elasticated, green, pack of 100



Bouffant Caps

Universal round cap made of light, air-permeable non
woven material.

- DH-12222* Blue, 50cm, pack of 100
DH-12214 Green, 56cm, pack of 100
DH-12404 White, 50cm, pack of 100



Patient Drapes

U-Drapes

3 layer


DH-301-628-S single (125 x 150cm)

DH-18-0101-10 pack of 10 mini split drapes for short procedures (58 x 80cm)

U-Drapes New Design

Two layer drape with redesigned aperture for greater comfort and ease of use.

DH-20-0366-S single (160 x 200cm)



Surgical Drapes

Non Adhesive

Two layer surgical drapes. High absorbent non woven laminated with PE-film on the bottom side. Drape is 100% fluid resistant. Sterile and packed separately.

DH-20-0131-S (40 x 60cm) pack of 10

DH-20-0103-S (75 x 75cm) pack of 10

DH-20-0105-S (75 x 120cm) pack of 10



Adhesive

Also available with adhesive edge for secure positioning

DH-20-0231-S (40 x 60cm) pack of 10

DH-20-0203-S (75 x 75cm) pack of 10

DH-20-0205-S (75 x 120cm) pack of 10



Tubing Sleeves

Standard Type

Full length tubes with adhesive strip

Pack of 20

DH-80-0103-S-20

Pack of 70

DH-80-0103-S

Concertina Type

Box of 15

DH-80-0125-S



Gloves Range

BODYGUARDS®
GLOVES RANGE

Handsafe Blue Nitrile Sterile Exam Gloves

- Contains no natural rubber latex proteins, easy donning
- Fit and comfort similar to natural latex
- More puncture resistance than latex
- Superior tactile sensitivity

GS-6905 small, 50prs

GS-6906 medium, 50prs

GS-6907 large, 50prs

LATEX FREE

TOP
SELLER!



Biogel Gloves

Sterile and hypoallergenic. Made from natural rubber latex. Powder-free surgical gloves with Biogel coating. Available in boxes of 10.

Sizes:

REG-002 6.0

REG-003 6.5

REG-004 7.0

REG-005 7.5

REG-006 8.0

POWDER FREE

Face Masks with loops

LHM-E1301 Box of 50 masks, Blue



Miscellaneous Items

Sterile Saline

BAXT-500

500ml saline semi-rigid pouch (box of 10)

BAXT-1000

1 litre plastic bottles (box of 8)

Sterile Water

BAXW-1000

1 litre plastic bottles (box of 8)

Spike Adaptor

BAXT-1000-S

Fits both saline and water bottles (box of 10)

Suction Tubing

PEN-001 (box of 20)

For connecting suction source to waste collection system, catheters, yankauer tips, probes and other devices. Sterile, double wrapped individual peel pouches. Latex free.

Mini-Yankauer Tips

PEN-002 (packs of 10)

Individually packaged clear sterile tube with vacuum control. Overall length 160mm.

Dispomedic Scrub

CV-31020 (box of 30)

Scrub brush with Chlorhexidine Gluconate 4% including nail cleaner.

Instant Cold Pack

DIS-001 (box of 25)

Disposable instant cold pack.

Squeeze to activate.

Comfort texture finish.

Convenient size 24 x12cm.

Steritips

Fine suction tips for bone collection and general use.

ROE-001 2.5mm steritip (blue tip)

20 per box including adaptor.

ROE-007 1.2mm steritip (green tip)

20 per box including adaptor.

Light Handle Covers

DH-25-0236-2-S Standard cover, in 2's, box of 100

Motor Cover

DH-71-0102-S 85 x 85cm elasticated, box of 15

DH-71-0501-S 70 x 80cm rectangular bag, box of 20

DH-71-0306-S 15 x 20cm Clearcover

DH-71-0707-5-10 27.5 x 20cm Clearcover

Sterile Towels

DH-98101-2-S box of 50 packs packed in 2's, 32 x 50cm

EZ Removable ink for Botox, dermal fillers and other non-surgical procedures



These regular tip markers have an easy to remove ink (use alcohol swab) and come in 3 colours

EZ-1440-30

Mini markers, white, regular tip, box of 30.

EZ-1445-30

Mini markers, red, regular tip, box of 30.

EZ-1444-30

Mini markers, green, regular tip, box of 30.

EZ-1440-240

Mini markers, white, regular tip, box of 240.

EZ-1445-240

Mini markers, red, regular tip, box of 240.

EZ-1444-240

Mini markers, green, regular tip, box of 240.

Swabs

Sterile gauze swabs in peelable blister packs (packed in 5's, all 8ply)

DH-91199-5-S Gauze swab 8ply 5 x 5cm (packed in 5's) x 300

DH-91200-5-S Gauze swab 8ply 7.5 x 7.5cm (packed in 5's) x 200

DH-91201-5-S Gauze swab 8ply 10 x 10cm (packed in 5's) x 200

Shoe Covers

Shoe cover - SPECIAL - length 36cms

DH-31501 pack of 100, Blue Shoe cover - SPECIAL - length 50cms

DH-30505 pack of 100, White

Over Mittens

DH-71-0123-5 pack of 50

Fluid Administration Sets

For most types of surgical units.

BIEN-IRRIGALINE Fluid set for the Bien Air Chiropro, ChiroproL and iChiropro (pack of 10)

BIEN-ICHIROLINES As above but dual end tube for internal/external irrigation handpieces (pack of 10)

N-1706 Fluid set for Nouvag models MD11/MD30 (pack of 10)

WH-04363600 Fluid set for ImplantMed (2.2m set, pack of 6)

Incise Foils

Clear dressing sheets for covering the patients face during surgery.

DH-76-0101-S 15 x 20cm (pack of 10)

Products	Page	Products	Page	Products	Page
Admin sets	83	Chisels	17	Instrument sets	36,37
Alveoplasty kit	14	Cold packs	83	Krekeler calipers	10
Assunyl sutures	63	Core tool	68	Light handle covers	82
Assupro sutures	62	Core tool - modified	69	Lindeman burs	66
Attachment processing material	72	Custom procedure packs	78	Locator overdenture implant system	66,67
Benex	20	Cytoplast dressings and membranes	56,57	Locator attachment systems	68-71
Bilumix headlamps and loupes	4,5	Cytoplast sutures	60	Locator Fixed system	69
Biogel gloves	81	Diamond discs and burs	14,16	Locator R-Tx	70,71
Biomaterial syringe	43	Drapes	80	Locator multi-unit abutment system	75
Bone block holders	11	Elevators	26,27	Loupes and headlamps	4,5
Bone mill	11	Explantation forceps	41	Luxators	22,23
Bone scrapers	34	Extraction pliers	24,25	Masks	81
Bone screw sets	13,18,19	Face masks	81	Mallets	17
Bone and socket sponges	44,56	Fast system plating kit	54	Master core system	15
Bone tack set	10,12,18	Fluid admin sets	83	Master Pin Control	12
Botox skin markers	83	Glicofil sutures	63	Membranes	46-59
Bowls and gallipots	39	Gloves	81	Membrane fixation system	10,12,18
Burs	14	Gowns - sterile	79	Micro screw system	13
Calipers	10	Haemostats	56	Microsaw scraper	34
Caps	79	Hats and head covers	79	Microsurgery instrument set	36
Cassettes	40	Headlamps and loupes	4,5	Mixing cup	39
Castroviejo needle holders	30	Hose covers	80	Monocryl sutures	61
Cement	73	Implant systems	64-67	Monofil sutures	63
Chairside attachment processing system	72	Improv implant cement	73	Monoglyc sutures	55
Choukroun PRF system	52	Incise foils	83	Motor covers	82

Products	Page	Products	Page	Products	Page
Multi-Unit abutment system	75	Roba extraction forceps	24,25	Terms and conditions	86
Needle Holders	30	Rongeurs	11	Tidespon sponge cubes	56
Nizam tissue punch	16	Root fragment forceps - Witzel	10	Tissue punch	15,16
O-Ring attachment system - Saturno	74	RPM reinforced membranes	58	Titanium mesh	56
Over mittens	83	RTM collagen membranes	56	Titanium reinforced membranes	58,59
Oxitamp haemostats	59	RTM wound dressings	56	Titanium trauma splints	36
Patient drapes	80	Safescrapers	34	Towels - sterile	82
Periosteal elevators	26,27	Saline	82	Trauma splints	36
Periotomes	21	SaturnO O-Ring attachment system	74	Trephines	15
Piezosurgery Touch	6	Scalpels and blades	33	Trolleys	6
Piezosurgery White	6	Scissors	28,29	Tubing	82
Piezosurgery tips and kits	7-9	Screws	13,19	Tubing sleeves	80
Plating system	54	Shoe covers	83	Tunnelling instruments	35
Prep and polish kit	72	Sinus lift sets	38,39	Tweezers	31
PRF system	52,53	Soft brushing kit	55	TXT-200 PTFE membranes	57
Procedure packs	76-78	Spike adaptors	82	Vicryl sutures	61
Profix bone screw system	18,19	Sterile packs	76-78	Water	82
Prolene sutures	61	Stiller sinus lift set	39	Wen tension release comb	55
Prosthetic recovery kit	37	Steritips	82	Witzel root fragment forceps	10
PTFE membranes	48-51,57-59	Suction tubing	82	Wound dressings	56
PTFE sutures	47,60	Sutures	55,	X-Tool luxators	22,23
Punches	16	Swabs	83	Yankauer tips	82
Purgo biomaterials and membranes	42-51	Tacks	10,12,18		
Retractors	32	Tension release comb	55		
REX Wedge implants	64,65	Tenting system	18,19		

1. (a) All contracts entered into between Jewel Management Ltd trading as General Medical ("the Seller") and any person firm or company ("the Buyer") purchasing goods and/or services from the Seller shall be subject to the following terms and conditions.
 (b) The Buyer's terms and conditions shall only be incorporated in the Contract:
 - (i) if they are not at variance with the Seller's conditions; and
 - (ii) if they are specifically incorporated into the Contract by the Seller in writing.
 If the Buyer's order purports to exclude these Conditions such clause shall be ineffective unless within three days of the date of the Contract (which shall be the date of the receipt by the Buyer of the Seller's written acceptance) the Buyer shall have obtained the Seller's express written agreement to vary these Conditions.
2. Unless otherwise stated all prices are net and are those ruling at the time of despatch. Published price lists and quoted prices while given in good faith are subject to change without notice due to suppliers' price fluctuations and other causes beyond the Seller's control. Quoted prices apply only to the stipulated quantities and do not necessarily hold for lesser quantities. All prices are subject where appropriate to Value Added Tax at the rate ruling at the time of despatch.
3. Payment is strictly net and due on delivery. The Seller reserves the right to charge interest on a day to day basis on accounts outstanding beyond the time specified in this Condition the rate of interest being 2 per cent per month until payment is received.
4. (a) Delivery and the cost of carriage shall be borne by the Buyer at their prevailing rates from time to time.
 (b) The Seller will use its best endeavours to meet quoted delivery dates but such dates are estimates only and the Seller shall not be liable for the consequences of any delay.
 (c) The Seller shall in any event be entitled to extend the delivery period where delay is caused by factors outside the Seller's reasonable control.
 (d) The Seller shall be entitled to deliver goods comprised in one order in several consignments in which event each such consignment shall be treated as the subject of a separate contract.
5. Damaged/faulty goods or short deliveries must be notified to the Seller within forty eight hours of delivery or where applicable installation and the goods and packaging material retained for inspection otherwise no liability can be accepted.
 The Seller's liability in respect of faulty goods shall be limited to giving the Buyer the benefit of any guarantee given by the manufacturer of such goods.
6. The title in the goods shall remain vested in the Seller and shall not pass to the Buyer until the Buyer has made payment in full of the purchase price. So long as the property in the goods remains vested in the Seller the Seller shall be at liberty at any time to retake possession thereof and for that purpose to enter upon any premises of the Buyer. Notwithstanding the foregoing the goods are at the entire risk of the Buyer from the point at which delivery is made.
7. Except where the goods are sold to a person dealing as a consumer (within the meaning of the Unfair Contract Terms Act 1977) all warranties conditions or other terms implied by statute or common law are excluded to the fullest extent permitted by law.
8. Orders placed cannot be cancelled except with the Seller's consent and on terms which will indemnify the seller against loss and expenses incurred. Any goods returned without the Seller's consent will not be accepted for credit. Where applicable goods will only be accepted for return within 30 days of the original date of invoice and will be subject to a minimum 10% restocking charge. The onus of proof of return remains with the customer. Under no circumstances are sterile products returnable.
9. Where stock is held specifically for a customer, all such stock will be deemed to be held on a consignment basis and can be delivered and invoiced at any time.
10. This contract shall be deemed to have been made in England the law applicable to it is the Law of England and Wales and the matter shall be tried in the Trowbridge County Court which will have exclusive jurisdiction in the case.
11. Finance options as shown in this catalogue are for illustrative purposes only and do not constitute a binding offer of finance. Finance facilities are provided by Medifinance Ltd who are authorised and regulated by the Financial Conduct Authority.



Order Online at www.generalmedical.co.uk
Or call 01380 734990 or email info@generalmedical.co.uk